

## LEXSYcon2.1 Expression Kit

**For constitutive expression of recombinant proteins in *Leishmania tarentolae***

**Cat. No. EGE-1310**

**contains one of the pLEXSY-2.1 vectors:**

**pLEXSY-ble2.1 (Cat. No. EGE-271)**

**pLEXSY-hyg2.1 (Cat. No. EGE-272)**

**pLEXSY-neo2.1 (Cat. No. EGE-273)**

**pLEXSY-sat2.1 (Cat. No. EGE-274)**

**pLEXSY-bsd2.1 (Cat. No. EGE-275)**

**pLEXSY-pac2.1 (Cat. No. EGE-276)**

**General purpose:**

**FOR GENERAL LABORATORY USE.**

**NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC**

**OR THERAPEUTIC USE.**

<b>1</b>	<b>INTRODUCTION</b>	<b>3</b>
<b>2</b>	<b>KIT COMPONENTS AND STORAGE CONDITIONS</b>	<b>4</b>
2.1	LEXSY host P10	4
2.2	pLEXSY-2.1 expression vector	4
2.3	Primers for diagnostic PCRs and sequencing	4
2.4	Ingredients for 1 litre of culture medium	5
2.5	Equipment and materials supplied by customer	5
<b>3</b>	<b>CULTURING OF LEXSY HOSTS AND EXPRESSION STRAINS</b>	<b>6</b>
3.1	Preparation of growth medium	6
3.2	Culturing conditions	6
3.3	Storage of LEXSY host and recombinants - Glycerol stocks preparation	7
<b>4</b>	<b>ENGINEERING OF LEXSY P10 EXPRESSION STRAINS</b>	<b>8</b>
4.1	Amplification of target gene	9
4.2	Insertion of target gene into pLEXSY-2.1 expression vector	11
4.3	Preparation of the expression plasmid	11
4.4	Transfection of the LEXSY host strain by electroporation	12
<b>5</b>	<b>SELECTION OF TRANSGENIC LEXSY STRAINS</b>	<b>13</b>
5.1	Clonal selection by plating on solidified media	13
5.2	Polyclonal selection in suspension culture	13
5.3	Confirmation of genomic integration by diagnostic PCR	14
5.4	Evaluation of target protein expression	15
<b>6</b>	<b>LICENSING INFORMATION</b>	<b>16</b>
<b>7</b>	<b>LITERATURE</b>	<b>16</b>
<b>8</b>	<b>APPENDIX</b>	<b>17</b>
8.1	Maps of the pLEXSY-2.1 expression vectors	17
8.2	Sequences of the primers available for LEXSYcon2 kits	23
8.3	Preparation of LEXSY BHI agar plates for clonal selection	26
8.4	Alternative electroporation protocols	27
8.5	How to grow a <i>Leishmania</i> culture	28
8.6	Correlation of optical density and cell concentration	29

# 1 INTRODUCTION


**LEXSY – the unique protein expression platform offered by Jena Bioscience** combines eukaryotic protein expression/folding/modification with robustness and easy handling. The unicellular protozoan host *Leishmania tarentolae* used in LEXSY was isolated from the Moorish gecko *Tarentola mauritanica* and kept in axenic culture over decades. It is not pathogenic to mammals and fully approved for use in biosafety level 1 (S1) laboratories.

## LEXSY features:

- Eukaryotic protein synthesis including chaperone system for correct folding (no inclusion bodies)
- Full range of Post-Translational Modifications including mammalian-type N-glycosylation, glypiation, phosphorylation, acetylation, prenylation, myristoylation, ADP-ribosylation, proteolytic processing, and oligomerization
- High expression-success rates with yields of up to 500 mg per litre of culture.

LEXSY is available in two principle configurations that are **constitutive** (Breitling *et al.* 2002) or **inducible** (Kushnir *et al.* 2005). In both configurations target proteins can be expressed either **intracellularly**, or be **secreted** into the culturing medium. An illustrated overview on LEXSY configurations available, features and applications can be downloaded from our website at [http://www.jenabioscience.com/images/b3e879b381/Lexsy\\_brochure\\_web.pdf](http://www.jenabioscience.com/images/b3e879b381/Lexsy_brochure_web.pdf).

The **LEXSYcon2.1 kit** contains one of the six expression vectors of choice which differ only in the antibiotic resistance marker: **pLEXSY-sat2.1** (encoding streptothricine acetyltransferase), **pLEXSY-hyg2.1** (encoding hygromycin phosphotransferase), **pLEXSY-ble2.1** (encoding bleomycin resistance gene), **pLEXSY-neo2.1** (encoding aminoglycoside phosphotransferase), **pLEXSY-bsd2.1** (encoding blasticidin-resistance gene) or **pLEXSY-pac2.1** (encoding puromycin N-acetyltransferase) allowing selection with the antibiotics **LEXSY NTC**, **LEXSY Hygro**, **LEXSY Bleo**, **LEXSY Neo**, **LEXSY BSD** or **LEXSY Puro** respectively. These expression vectors are designed for integration into the chromosomal 18SrRNA (*ssu*) locus of *L. tarentolae* (Figure 1). From this locus the target genes are transcribed by the strong RNA polymerase I of the host cells under the control of a chromosomal ribosomal promoter. Due to the positioning of an **additional *Leishmania* ribosomal promoter** in front of the expression cassette on the **new pLEXSY-2.1 vector family** the protein production potential of the previous pLEXSY-2 vectors could be significantly enhanced by 50%.



**Figure 1:** Architecture of the constitutive integrative LEXSY expression system. Target genes are inserted into the expression site of the integrative pLEXSYcon2.1 expression vector and, following linearisation of the obtained expression plasmid, introduced into the *Leishmania tarentolae* expression host. Ssu is the chromosomal 18S rRNA integration locus for the expression constructs. 5'ssu and 3'ssu are the ssu regions of the expression vector for homologous recombination with the host chromosome. Pr=Leishmania ribosomal promoter on the host chromosome and in front of the vector expression cassette.

## 2 KIT COMPONENTS AND STORAGE CONDITIONS

The kit is shipped on dry ice. Upon arrival the kit components should be stored at the appropriate temperature as indicated below.

### 2.1 LEXSY host P10

The kit contains 3 vials with 1.6 ml each of frozen glycerol stocks of **LEXSY host P10**. These stocks can be stored at -80°C for at least 1 year. However, we have successfully reactivated frozen LEXSY stocks after >10 years of storage at -80°C. For reactivation protocol refer to chapter 3.4.

### 2.2 pLEXSY-2.1 expression vector

*one of the vectors*

- **pLEXSY-ble2.1** (Cat. No. EGE-271) or
- **pLEXSY-hyg2.1** (Cat. No. EGE-272) or
- **pLEXSY-neo2.1** (Cat. No. EGE-273) or
- **pLEXSY-sat2.1** (Cat. No. EGE-274) or
- **pLEXSY-bsd2.1** (Cat. No. EGE-275) or
- **pLEXSY-pac2.1** (Cat. No. EGE-276)
  - 5 µg in 10 mM TrisHCl pH 8.0
  - store at -20°C

For vector maps see Appendix 8.1. The DNA sequences can be downloaded from the website of JBS at [http://www.jenabioscience.com/cms/en/1/catalog/1856\\_vectors\\_and\\_primers.html](http://www.jenabioscience.com/cms/en/1/catalog/1856_vectors_and_primers.html).

### 2.3 Primers for diagnostic PCRs and sequencing

*one of the vector-specific primers*

- **ble** forward primer **A708**, for **pLEXSY-ble2.1** or
- **hyg** forward primer **A3804**, for **pLEXSY-hyg2.1** or
- **neo** forward primer **H4662**, for **pLEXSY-neo2.1** or
- **sat** forward primer **F2999**, for **pLEXSY-sat2.1** or
- **bsd** forward primer **H4895**, for **pLEXSY-bsd2.1** or
- **pac** forward primer **H4896**, for **pLEXSY-pac2.1**

*and the primers*

- Insert sequencing forward primer **P1442**
- Insert sequencing reverse primer **A264**
- 5`utr (*aprt*) reverse primer **A1715**
- **ssu** forward primer **F3001**
- **ssu** reverse primer **F3002**
  - 50 µl each at 50 µM in 10 mM TrisHCl pH 8
  - store at -20°C

For primer sequences see Appendix 8.2.

## 2.4 Ingredients for 1 litre of culturing medium

**LEXSY BHI**, Powder for preparation of liquid culturing medium (Cat. No. ML-412S)

- 37 g (2 x 18.5 g for 500 ml each)
- store at ambient temperature
- stable for 12 months.

### Additives

- **Hemin** 0.25% porcine Hemin in 30% Triethanolamine (Cat. No. ML-108S)
  - 2 ml, ready-to-use 500x stock solution, filter sterilized
  - store at 4°C in the dark (Hemin is light sensitive)
  - stable for 12 months.
- **Pen-Strep** 10.000 units of penicillin (base) and 10.000 µg of streptomycin (base)/ml as penicillin G sodium salt and streptomycin sulfate in 0.85% saline (Cat. No. ML-105S) may be added to avoid bacterial contaminations
  - 5 ml, ready-to-use 200x stock solution, filter sterilized
  - store at -20°C
  - stable for 12 months.

### one of the selection antibiotics

**LEXSY Bleo** (Cat. No. AB-103S) for **pLEXSY-ble2.1**, ready-to-use 1000x stock solution, 100 mg/ml or

**LEXSY Hygro** (Cat. No. AB-104S) for **pLEXSY-hyg2.1**, ready-to-use 1000x stock solution, 100 mg/ml or

**LEXSY Neo\*** (Cat. No. AB-105S) for **pLEXSY-neo2.1**, ready-to-use 1000x stock solution, 50 mg/ml or

**LEXSY NTC** (Cat. No. AB-101S) for **pLEXSY-sat2.1**, ready-to-use 1000x stock solution, 100 mg/ml or

**LEXSY BSD** for **pLEXSY-bsd2.1**, ready-to-use 1000x stock solution, 10 mg/ml or

**LEXSY Pac** for **pLEXSY-pac2.1**, ready-to-use 1000x stock solution, 10 mg/ml

- 1 ml each, filter sterilized
- store at -20°C
- stable for 12 months.

\* Please note, that **LEXSY Neo** (JBS Cat. No. AB-105) is **not the antibiotic Neomycin** which must not be used for selection in LEXSY.

## 2.5 Equipment and materials supplied by customer

- Incubator at 26°C (no CO<sub>2</sub> incubator required), inverse (or standard) microscope
- Electroporation device, e.g. BioRad GENEPULSER II with PULSE CONTROLLER II and CAPACITANCE EXTENDER PLUS or GENE PULSER Xcell with PC and CE Modules, or Multiporator, or Nucleofector
- Cooling and freezing capacities at +4°C, -20°C and -80°C
- Standard molecular biology equipment for PCR, cloning and protein analysis.

### 3 CULTURING OF LEXSY HOST P10 AND EXPRESSION STRAINS

The standard culturing medium for strain maintenance, transfection, cryoconservation and expression evaluation is **LEXSY BHI Medium** (Cat. No. ML-411, 412). It must be supplemented with **Hemin**, which is essential for *Leishmania*. There is no need to add sera to complex media. Addition of fetal calf serum will not enhance growth of *L. tarentolae* in complex media. To prevent bacterial infections, Penicillin and Streptomycin (**Pen-Strep**) may be added.

#### 3.1 Preparation of LEXSY BHI Medium

Dissolve 37 g/l **LEXSY BHI** powder in de-ionized water and autoclave **exactly** 15 min at 121°C. Control temperature profile with a thermosensor in a reference flask with the same volume liquid. Note that overexposure may result in decomposition of media ingredients (e.g. indicated by a dark brown color) which adversely affects growth of LEXSY strains. Properly autoclaved medium is of amber color. Store core medium at room temperature before addition of Hemin and PenStrep (stable for 12 months).

Add to 500 ml **LEXSY BHI Medium**

- 2.5 ml of 200x **Pen-Strep** stock solution
- 1.0 ml of 500x **Hemin** stock solution (final concentration 5 µg/ml)
- store at 4°C in the dark and use within 2 weeks after supplementation. If the completed media are to be used after this period, appropriate amounts of additives have to be re-added.

#### 3.2 Culturing conditions

All culturing is carried out **at 26°C** in the **dark** under **aerated** conditions (no CO<sub>2</sub> incubator required). **Hemin** must be added to all culturing media.

**Standard culturing for strain maintenance** is performed as continuous static suspension cultures in ventilated tissue culture (TC) flasks positioned upright with regular successive dilutions into **LEXSY BHI Medium**. Best results are obtained with dilutions at early stationary phase (OD 1.4-2 corresponding to ca. 6-8x 10<sup>7</sup> cells/ml under these conditions). We find it convenient to dilute 10 ml cultures 1:50 on Monday and 1:20 on Friday. Avoid repeated successive dilution of cultures of lower cell densities as this may reduce growth. However, occasional higher dilutions of stationary cells at e.g. 1:100 will not adversely affect subsequent growth. Do not use agitation for strain maintenance (ref. to Appendix 8.5).

To prevent genomic changes upon long-term culturing of P10 host it is recommended, not to exceed successive passages for more than three months. Instead, a **new culture of the P10 host should be started from a glycerol stock every three months**.

For **culturing for transfection** refer to chapter 4.4 and for **protein expression** see chapter 5.3

To conveniently monitor the growth of the LEXSY cultures calibrate the OD readings of your spectrophotometer at a defined wavelength between 550 and 600 nm to the cell densities at different time points during growth of a suspension culture by taking OD readings and counting cells from the same sample e.g. in a particle counter or in a Neubauer chamber after immobilisation in 3% (final concentration) formalin. Dependent on the spectrophotometer and wavelength used, this correlation may be different from the data specified in this manual (Figure 12 in Appendix 8.6).

In case you encounter growth problems with the host or LEXSY expression strains, centrifuge cells 3 min at 2000 x g, resuspend pellet carefully in fresh medium and continue incubation.

Culturing may be performed in

- ventilated tissue culture (TC) flasks for suspension cultures, culture volume 10-100 ml
- Erlenmeyer flasks, agitated at approx. 120-140 rpm, culture volume 50-200 ml
- Fernbach flasks, agitated at approx. 90 rpm, culture volume ca. 0.1-0.5 litre
- standard bioreactors, culture volume 1-100 litres at 100-200 rpm and  $pO_2 > 40\%$ .

### 3.3 Storage of LEXSY host and recombinant strains by cryoconservation

The LEXSY host and LEXSY expression strains may be stored at  $-80^\circ\text{C}$  in 10-20% glycerol for at least 12 months. However, we recovered viable cells under these conditions after >10 years of storage, without any loss of vitality. Glycerol stocks may be prepared from 1:10 inoculated static cultures after 24 h if TC flasks are incubated flat at  $26^\circ\text{C}$  (OD at 600 nm ca. 1.4 corresponding to ca.  $5 \times 10^7$  cells/ml). Prior to conservation, check vitality of cells by microscopy. Cells should be motile and drop-like (Figure 11, Appendix 8.5) but not of needle-like appearance. Do not use other media than **LEXSY BHI Medium** for cryoconservation.

#### Glycerol stocks preparation using a cryocontainer

Using this method, freezing proceeds continuously with  $-1^\circ\text{C}/\text{min}$

- Add 1.2 ml autoclaved Glycerol (80% by weight per volume) to a sterile 15 ml Falcon tube
- Withdraw 3.6 ml of culture OD ca. 1.4-2 (ca.  $5 \times 10^7$  cells/ml)
- Mix with glycerol and distribute 3 x 1.6 ml each to sterile cryovials
- keep 10 min at room temperature
- transfer to a cryocontainer at  $4^\circ\text{C}$  containing fresh isopropanol
- keep 10 min at  $4^\circ\text{C}$
- transfer to  $-80^\circ\text{C}$  over night
- distribute to storage box for long term storage.

#### Glycerol stocks preparation by stepwise cool-down protocol

Using this method, freezing proceeds in steps of  $0^\circ\text{C}/-20^\circ\text{C}/-80^\circ\text{C}$

- Add 1.2 ml autoclaved Glycerol (80% by weight per volume) to a sterile 15 ml Falcon tube
- Withdraw 3.6 ml of culture OD ca. 1.4-2 (ca.  $5 \times 10^7$  cells/ml)
- Mix with glycerol and distribute 3 x 1.6 ml each to sterile cryovials
- keep 10 min at room temperature
- keep 1 h on wet ice
- keep o/n at  $-20^\circ\text{C}$
- transfer to  $-80^\circ\text{C}$  for long term storage.

Both protocols are tolerated well by *L. tarentolae*. However, to avoid loss of strains it is recommended to check the reactivation of one sample of the prepared batch of glycerol stocks prior to stopping respective suspension culturing of parent culture.

#### Glycerol stocks reactivation

- Thaw frozen glycerol stock on ice (ca. 20 min)
- Inoculate the **entire content** of the vial into 10 ml of **LEXSY BHI Medium** with appropriate antibiotic(s). Motile cells can be observed immediately after inoculation by microscopy
- Incubate as static suspension culture in ventilated TC flask (flat) dark at 26°C until culture gets turbid (OD 1.4-2; ca. 6-8 x10<sup>7</sup> cells/ml). This usually takes 2 days; wait longer if cells recover more slowly and follow status by microscopy
- Dilute dense culture 1:10-1:20 into fresh LEXSY BHI Medium and incubate for 3 days. Do not dilute culture of low density. For strain maintenance dilute into fresh LEXSY BHI Medium on Monday and Friday each week (see chapter 3.2).

#### 4 ENGINEERING OF LEXSY P10 EXPRESSION STRAINS

This kit was designed for constitutive expression of target proteins following integration of the expression cassette into the chromosomal 18S rRNA locus (*ssu*) making use of the high level transcription rates by the host RNA polymerase I (Figure 1).


Due to the positioning of an **additional *Leishmania* ribosomal promoter** in front of the expression cassette on the **new pLEXSY-2.1 vector family** the protein production potential of the previous pLEXSY-2 vectors could be significantly enhanced by 50%.

In the first step the ORF for the target gene is supplied with linker sequences containing restriction sites allowing insertion into the **pLEXSY-2.1 vectors** (Figure 2). These vectors contain optimized non-translated regions flanking the target gene insertion sites, which provide the splicing signals for posttranscriptional mRNA processing. Dependent on the cloning strategy, target proteins are destined either for intracellular or for secretory expression (Figure 2). Secretory expression is achieved by fusion of the target ORF to the *Leishmania* signal peptide coding region present on the vector. For cloning traditional ligase-dependent or alternative ligase-free methods may be followed. The ligase-free **CyClone** technology of Jena Bioscience (JBS Cat. No. PP-110) offers the advantage of online design and ordering of primer pairs for insert amplification with vector-compatible overhangs using the free Primer Designer software of Jena Bioscience [http://www.jenabioscience.com/cms/en/1/catalog/1856\\_vectors\\_and\\_primers.html](http://www.jenabioscience.com/cms/en/1/catalog/1856_vectors_and_primers.html).

Following plasmid construction in *E. coli* the LEXSY host P10 is transfected with the linearized expression construct by electroporation. Clonal or polyclonal selection strategies can be applied for establishment of recombinant strains.







**Figure 2:** Map of the **pLEXSY-2.1** vector family with cloning sites for the target genes replacing the 1 kb stuffer fragment. **5'ssu** and **3'ssu** are regions for homologous recombination into the host chromosome following linearization of the expression plasmid with **SwaI**. **Pr** is a *Leishmania* ribosomal promoter in front of the vector expression cassette. **Utr1** derived from 0.4k-IR of *L. tarentolae aprt*, **utr2** from 1.4k-IR *camCB* and **utr3** from 1.7k-IR are optimized non-translated gene-flanking regions providing the splicing signals for posttranscriptional mRNA processing for expression of target and marker genes in the LEXSY host P10. **SP** designates the signal peptide of *L. mexicana* secreted acid phosphatase LMSAP1 (Wiese *et al.* 1995), **SPCS** the signal peptide cleavage site and **H6** the hexa-Histidine stretch (underlined in upper sequence). Alternative cloning strategies result in cytosolic (**c**) or secretory (**s**) expression of the target protein. Please, note that the *KasI* site at the SPCS is not unique and not suited for target gene insertion. As markers are available the **ble** (bleomycin resistance), **hyg** (hygromycin phosphotransferase), **neo** (aminogluco-side phosphotransferase) **sat** (streptothricine acetyltransferase), **bsd** (blastidicin-resistance) or **pac** (puromycin N-acetyltransferase) genes. The DNA sequences of all pLEXSY-2.1 vectors can be downloaded from the website of JBS at [http://www.jenabioscience.com/cms/en/1/catalog/1856\\_vectors\\_and\\_primers.html](http://www.jenabioscience.com/cms/en/1/catalog/1856_vectors_and_primers.html).

## 4.2 Insertion of target gene into pLEXSY-2 expression vector

The **pLEXSY-2 vectors** allow directional insertion of the target gene cassette by replacement of a 1 kbp stuffer fragment (Figure 2). This is advantageous to insertion into a multiple cloning site since vector cleavage by both restriction enzymes can be monitored by the appearance of the stuffer fragment.

- Digest the expression vector provided in the kit with the appropriate restriction enzymes.
- Gel-isolate the large 7-8 kbp fragment and prepare it for ligation.
- Ligate vector and target gene from chapter 4.1 with T4 DNA ligase (JBS Cat. No. EN-149). For discrimination against the original expression vector you may treat the ligation mix with HpaI after ligation reaction if your target gene does not contain this site. Alternatively, perform **CyClone** reaction if ligase-free method was chosen (JBS Cat. No. PP-110)
- Transform competent *E. coli* cells which tolerate Leishmania sequences as XL-10, DH10B, Stbl2, Stbl4, SURE, DH5α etc.
- Select recombinant *E. coli* clones with ampicillin at 30°C and screen for the presence of the insert in the plasmids. Insert screening may be performed by colony PCR or restriction analysis of recombinant plasmids isolated from a small number of cultures grown o/n in 1 - 3 ml LB with ampicillin at 30°C. We recommend incubation of all *E. coli* strains with pLEXSY plasmids **at 30°C and not at 37°C** for plasmid stability reasons
- Prepare at least 10 µg plasmid DNA from a positive clone for sequence confirmation, restriction and subsequent transfection. Usually, it is sufficient to isolate plasmid DNA with a commercial kit from 50 ml o/n culture grown at 30°C.
- Confirm plasmid identity by restriction and sequence analysis of the insert and of the vector/insert fusions using the forward **P1442** and reverse **A264** sequencing primers included in the kit. Both primers are proved for function in cycle sequencing protocols. Forward primer **P1442** anneals 5' of the insert (ca. 80 bp in front of ATG), reverse primer **A264** 3' of the insert (ca. 80 bp after the stop codon). The primer sequences are shown in Appendix 8.2 of this manual, the DNA sequence of the vector including description can be downloaded as .txt file from the LEXSY section of our website [http://www.jenabioscience.com/cms/en/1/catalog/1862\\_lexsy\\_expression\\_vectors.html](http://www.jenabioscience.com/cms/en/1/catalog/1862_lexsy_expression_vectors.html).

## 4.3 Preparation of the expression plasmid for LEXSY host transfection

- Digest to completion with **SwaI** ca. 10 µg of the obtained expression plasmid containing the target gene. This treatment will generate a 2.9 kbp fragment representing the *E. coli* part and a larger fragment representing the linearized expression cassette with the target gene to be integrated into the chromosomal *ssu* locus of the LEXSY host
- For best performance, gel-isolation of the expression cassette with an Agarose Gel Extraction Kit (Cat. No. PP-202) is recommended (but optional)
- If the digested plasmid is used without fragment purification, enzymes and buffer salts may be removed with a PCR Purification Kit (Cat. No. PP-201). Alternatively, precipitate the digest with ethanol, wash with 70% ethanol and re-dissolve in max. 50 µl sterile double distilled water or 10 mM Tris pH 8 per transfection. Control the quality of the digest and the DNA concentration by gel electrophoresis. This preparation is now ready for transfection.

#### 4.4 Transfection of the LEXSY host strain P10 by electroporation

For efficient transfection it is recommended to prime the LEXSY host by successively transferring the cells to fresh medium at 1:20 to 1:50 dilutions twice a week (refer to chapter 3.2). Do not use the first inoculation culture from a glycerol stock immediately for transfection, but passage the culture several times before electroporation as described above. However, use for transfection a culture which has been passaged this way for less than three months (refer to chapter 3.2).

- On Friday inoculate *L. tarentolae* pre-culture 1:20 in 10 ml **LEXSY BHI** medium supplemented with Hemin and PenStrep (refer to chapter 3.1) and incubate in tissue culture (TC) flask upright @ 26°C dark until Monday
- On Monday dilute pre-culture 1:10 into 10 ml of same medium and incubate in TC flask flat @ 26°C o/n
- On Tuesday check cell density of the culture until approx.  $6 \times 10^7$  cells/ml are reached (OD 1.4)\* and ensure by microscopy that the cells are vital and of droplike shape (Figure 11, Appendix 8.5)
- Spin cells 3 min, 2000g at room temperature and remove ½ volume of supernatant
- Resuspend pellet in remaining medium to get  $10^8$  cells/ml and put on wet ice for 10 min
- Have ready on wet ice in parallel tubes with 1-10 µg transforming DNA in max. 50 µl of water or 10 mM Tris buffer pH 8.0 and electroporation cuvettes d=2 mm\*\*
- Add 350 µl pre-chilled cells to the tube with DNA and transfer all 400 µl to the electroporation cuvette on wet ice
- Electroporate @ 450V, 450 µF and monitor pulse time (ca. 5-6 msec)\*\*\*
- Put cuvette back on ice for exactly 10 min
- Transfer electroporated cells with a capillary to 10 ml LEXSY BHI Medium supplemented with Hemin and PenStrep in a ventilated TC flask
- Incubate o/n @ 26°C flat as static suspension culture (ca. 20h to OD 0.3-0.4)
- Proceed with clonal or polyclonal selection (chapters 5.1 and 5.2).

\* if the cell density differs from this value, concentrate cells in the next step to get  $10^8$  cells/ml. For transfection, cultures between OD 1.0-1.8 can be used. Do not transfect cells if they are long and thin by microscopy.

\*\* use electroporation cuvettes with long electrodes. The entire volume of 0.4 ml must be between the electrodes. Do not use electroporation cuvettes with short electrodes leaving most of the volume outside of the linear electric field.

\*\*\* using BioRad GENEPULSER II with PULSE CONTROLLER II and CAPACITANCE EXTENDER PLUS or GENE PULSER Xcell with PC and CE Modules. The resistance of the sample is 20 Ohms. With GENE PULSER Xcell you may alternatively use the **Time constant protocol** with the **settings 450 V and 3.5 ms** (Figure 10 in Appendix 8.4). Appendix 8.4 describes also an alternative **High voltage protocol** for transfection of LEXSY.

The protocol for electroporation of LEXSY with a **Multiporator** can be downloaded from the website of JBS at [http://www.jenabioscience.com/cms/en/1/browse/1888\\_lexsy\\_methods.html](http://www.jenabioscience.com/cms/en/1/browse/1888_lexsy_methods.html).

The protocol for electroporation of LEXSY with a **Nucleofector** is described in Vainio *et al.* 2009.

## 5 SELECTION OF TRANSGENIC LEXSY STRAINS

For establishing expression strains you have the option of clonal and/or polyclonal selection. We repeatedly found similar expression levels when comparing cultures derived from clonal (chapter 5.1) or polyclonal (chapter 5.2) selection strategies following transfection with linearized expression cassettes designed for chromosomal integration.

### 5.1 Clonal selection by plating on solidified media

For plating, LEXSY BHI agar plates are always freshly prepared on the day of plating as described in Appendix 8.3. For customer convenience, all components required for preparing LEXSY BHI agar plates are included in our **LEXSY Plating Kit** (Cat. No. ML-451).

- Pellet the cells from the transfected 10 ml o/n culture (chapter 4.4) for 3 min at 2000g and room temperature
- Remove the supernatant and resuspend the cells in approx. 50-100 µl of residual medium
- Carefully spread the resuspended cells on top of the freshly prepared selective LEXSY BHI agar
- Seal plates with parafilm and incubate bottom up at 26° for 5-7 days until small, defined colonies begin to appear (approx. 5-7 days after plating)
- Let them grow up to 1–2 mm diameter (approx. 7-9 days after plating)
- Transfer the selected colonies with a pipette tip into 1 ml of medium supplemented with the selective antibiotic in 24 well format and incubate as static or agitated\* suspension cultures. This should be done within three weeks post plating
- Proceed to 5.3 for confirmation of genomic integration.

\* Agitated suspension culturing in 24 well format may be performed on a microplate shaker with high amplitude (e.g. DESAGA TPM-2).

### 5.2 Polyclonal selection in suspension culture

- As soon as the 10 ml o/n cultures obtained from the transfection experiments (chapter 4.4) start to get slightly turbid ( $OD_{600}$  0.4; ca.  $10^7$  cells/ml; usually approx. 20 h after electroporation), add the appropriate selective LEXSY antibiotic to the recommended final concentration from the filter-sterilized stocks provided in the kit (chapter 2.4) and continue incubation as static suspension culture at 26°C. **Don't let the cultures overgrow before selection, since it will take longer to kill non-recombinant cells!**
- Follow the status of the cultures microscopically and visually until you start seeing the difference to the cells electroporated without DNA under the same conditions (negative control). Recombinant cells are motile, of drop-like shape and grow as a "cloudy" suspension culture whereas the cells in the negative control begin to die during the selection period and appear as spherical or irregular forms without flagella under the microscope. Visually, however, the negative control may also appear as a turbid suspension
- It usually takes one consecutive transfer into fresh medium with the appropriate antibiotics at 1:10 inoculation rate to get a turbid culture of the antibiotic-resistant recombinant cell line and a clear negative control (absence of growth). This passage should be performed within 7 days after first drug addition, usually at day 5. Do not wait for a longer period to passage the culture even if the negative control also got turbid
- if you feel that the transfected cultures do not start to grow under selection within the first week post electroporation, centrifuge cells 3 min at 2000g, resuspend pellet carefully in fresh medium

with selective antibiotics and continue incubation. Also, you may combine 1:10 dilution with medium renewal of the remaining culture at day 5 of selection

- If at the time of first addition of selective antibiotic the cultures were already too dense, the negative control may not be dead after the first 1:10 passage. In this case another (or even more) subsequent 1:10 dilution(s) into fresh medium with selective antibiotic may be needed to get a turbid recombinant culture and a clear negative control

An illustration of the selection procedure can be downloaded from our website at [http://www.jenabioscience.com/images/b3e879b381/illustration\\_LEXSY\\_Transfection\\_Selection.pdf](http://www.jenabioscience.com/images/b3e879b381/illustration_LEXSY_Transfection_Selection.pdf).

### 5.3 Confirmation of genomic integration by diagnostic PCR

Integration of the expression cassette into the *ssu* locus can be confirmed by diagnostic PCR using genomic DNA of transgenic strains as template. For this purpose primer pairs including one primer hybridizing within the expression cassette and one primer hybridizing to a chromosomal *ssu*-flanking sequence not present on the plasmid are used (Figure 9 of Appendix 8.2). A set of primers for such diagnostic PCRs as well as for insert sequencing are included in each kit. The primer sequences are shown in appendix 8.2 of this manual.

- Prepare genomic DNA from 2 ml of a dense culture (OD approx. 2-3) by conventional phenol/chloroform extraction or with a commercial kit (e.g. Cat-No. PP-206, 208). There is no need to add Lysozyme or Proteinase K during purification
- Perform diagnostic PCR with the primer pairs and under the conditions outlined in the table below

Control region	Forward primer	Reverse primer	Fragment size	Annealing temp.
<b><i>5'ssu - utr1</i></b>	F3001	A1715	1.6 kbp	60°C
<b><i>ble - 3'ssu</i></b>	A708	F3002	2.2 kbp	54°C
<b><i>hyg - 3'ssu</i></b>	A3804	F3002	1.8 kbp	53°C
<b><i>neo - 3'ssu</i></b>	H4662	F3002	2.6 kbp	53°C
<b><i>sat-3'ssu</i></b>	F2999	F3002	2.3 kbp	53°C
<b><i>bsd-3'ssu</i></b>	H4895	F3002	2.2 kbp	53°C
<b><i>pac-3'ssu</i></b>	H4896	F3002	2.4 kbp	53°C

- The PCR reactions will result in a characteristic fragment for each vector (see table above and Figure 9 of Appendix 8.2), which is not observed in control reactions where the template is the expression plasmid or genomic DNA from the LEXSY host. Additional diagnostic PCR reactions including target gene-specific primers may be performed.

### 5.4 Evaluation of target protein expression

Expression of the target protein in recombinant LEXSY strains may be evaluated by SDS-polyacrylamide gel electrophoresis (SDS-PAGE) and Western blotting of cell extracts or, in case of secretory expression, aliquots from supernatants. For obtaining optimal expression we recommend to check different culturing conditions and time of harvest for each individual protein.

Maintain selected recombinant clones from 5.1/5.2 as continuous 10 ml suspension cultures as described in chapter 3.2. Always add selective LEXSY antibiotics to these cultures. Antibiotics may be omitted from larger cultures destined for downstream processing.

- Inoculate 10-100 ml of LEXSY BHI medium 1:10 with the expression clone and grow at 26°C in TC flask (flat) as static suspension culture or as agitated culture on a microplate shaker with high amplitude (e.g. DESAGA TPM-2). Alternatively, grow culture in Erlenmeyer flasks at approx. 140 rpm. For culturing of larger volumes use Fernbach flasks (200-700 ml) or fermenters (>1 litre). For higher cell densities and yields, add 2 g/L glucose at OD>2
- Take aliquots at different time points of culturing (e.g. OD 2 - 10) for estimation of optimal harvest time. For the proteins tested we found harvest of agitated batch cultures ca. 48h post inoculation and of static cultures ca. 72h post inoculation were optimal
- For **analysis of intracellular expression** calculate volume of aliquot  $V \text{ [ml]} = 2/OD$ , e.g. withdraw 4 ml culture @ OD 0.5 or 2 ml culture @ OD 1.0 or 1 ml culture @ OD 2.0 etc.
  - Sediment cells 5 min at 3.000g
  - Resuspend pellet in 0.2 ml of gel loading buffer and apply 10-20 µl per lane on SDS-PAGE for Coomassie staining and/or Western blotting
- For **analysis of secretory protein expression** concentrate culture supernatants 100x with trichloroacetic acid (TCA) as follows:
  - Sediment cells from 10 ml of culture 10 min @ 3000g
  - Add 8 ml of (steril-filtered) supernatant to 2 ml of 50% ice-cold TCA to a final concentration of 10%. Sterilfiltration of supernatants prior to TCA precipitation avoids carry-over of cells and is optionally
  - Leave on ice for 30 min, spin 15 min 15.000g at 4°C
  - Remove supernatant completely and collect for safe waste disposal
  - Resuspend pellet in 1 ml of 80% acetone and transfer to an Eppendorf tube (the acetone-wash is performed to remove residual TCA)
  - Spin 15 min 15.000g at 4°C, aspirate supernatant and resuspend pellet in a final volume of 80 µl gel loading buffer (corresponding to 100x concentration)
  - Apply 10-20 µl sample/slot for SDS-PAGE and Western blotting
- For purification of target proteins use appropriate culturing in larger volumes as described in chapter 3.2 and use affinity or conventional techniques.

## 6 LICENSING INFORMATION

Purchase of the **LEXSY Expression Kits** includes a non-exclusive and non-transferable license for non-commercial research. Commercial use of the **LEXSY expression system**, however, requires separate licensing. Commercial use includes but is not limited to:

- the use of any protein or other substance produced by LEXSY as reagents in screening to discover and/or promote candidate compounds for sale to a customer, distributor, wholesaler or other end user in therapeutic, diagnostic, prophylactic, and/or veterinary areas
- the manufacture, sale or offer to sell of any product containing proteins or other substances produced by LEXSY
- the large-scale production of recombinant protein pharmaceuticals
- "Contract research" to any third party or "Contract manufacturing" for any third party that has not been granted a license to use LEXSY

Please, contact us at [expression@jenabioscience.com](mailto:expression@jenabioscience.com).

## 7 LITERATURE

Breitling R, Klingner S, Callewaert N, Pietrucha R, Geyer A, Ehrlich G, Hartung R, Müller A, Contreras R, Beverley S and Alexandrov K (2002) Non-pathogenic trypanosomatid protozoa as a platform for protein research and production. *Protein Expression and Purification* **25**: 209-218

Klatt S and Konthur Z (2012) Secretory signal peptide modification for optimized antibody-fragment expression-secretion in *Leishmania tarentolae*. *Microbial Cell Factories* **11**: 97

Kushnir S, Gase K, Breitling R and Alexandrov K (2005) Development of an inducible protein expression system based on the protozoan host *Leishmania tarentolae*. *Protein Expression and Purification* **42**: 37-46

Lukeš J, Paris Z, Regmi S, Breitling R, Mureev S, Kushnir S, Pyatkov K, Jirků M, and Alexandrov K (2006) Translational initiation in *Leishmania tarentolae* and *Phytomonas serpens* (Kinetoplastida) is strongly influenced by pre-ATG triplet and its 5' sequence context. *Molecular & Biochemical Parasitology* **148**: 125-132

Robinson KA and Beverley SM (2003) Improvements in transfection efficiency and tests of RNA interference (RNAi) approaches in the protozoan parasite *Leishmania*. *Molecular & Biochemical Parasitology* **128**: 217-228

Vainio S, Genest PA, Ter Ried B, van Luenen H and Borst P (2009) Evidence that J-binding protein 2 is a thymidine hydroxylase catalyzing the first step in the biosynthesis of DNA base J. *Molecular & Biochemical Parasitology* **164**: 157-161


Wiese M, Ilg T, Lottspeich F and Overath P (1995) Ser/Thr-rich repetitive motifs as targets for phosphoglycan modifications in *Leishmania mexicana* secreted acid phosphatase. *EMBO Journal* **14**: 1067-10748.

Please, visit [http://www.jenabioscience.com/cms/en/1/browse/1853\\_scientific\\_literature.html](http://www.jenabioscience.com/cms/en/1/browse/1853_scientific_literature.html) for further literature on LEXSY.




## 8 APPENDIX


### 8.1 Maps of the pLEXSY-2.1 expression vectors




**Figure 3:** Map of **pLEXSY-ble2.1** (Cat. No. EGE-271) expression vector with cloning sites for the target genes replacing the 1 kb stuffer fragment. **5'ssu** and **3'ssu** are regions for homologous recombination into the host chromosome following linearization of the expression plasmid with **Swal**. **Pr** is a *Leishmania* ribosomal promoter in front of the vector expression cassette. **Utr1** derived from 0.4k-IR of *L. tarentolae aprt*, **utr2** from 1.4k-IR *camCB* and **utr3** from 1.7k-IR are optimized gene-flanking non-translated regions providing the splicing signals for posttranscriptional mRNA processing for expression of target and marker genes in the LEXSY host P10. **SP** designates the signal peptide of *L. mexicana* secreted acid phosphatase LMSAP1 and **H6** the hexa-Histidine stretch. The sequence of the multiple cloning sites is indicated in Figure 2. The DNA Sequence and further description of pLEXSY-ble2.1 can be downloaded from our website at <https://www.jenabioscience.com/images/103bb272b3/EGE-271.txt>. Please note, that the *KasI* site at the SPCS is not unique and not suited for target gene insertion.




**Figure 4:** Map of **pLEXSY-hyg2.1** (Cat. No. EGE-272) expression vector with cloning sites for the target genes replacing the 1 kb stuffer fragment. **5' ssu** and **3' ssu** are regions for homologous recombination into the host chromosome following linearization of the expression plasmid with **Swal**. **Pr** is a *Leishmania* ribosomal promoter in front of the vector expression cassette. **Utr1** derived from 0.4k-IR of *L. tarentolae aprt*, **utr2** from 1.4k-IR *camCB* and **utr3** from 1.7k-IR are optimized gene-flanking non-translated regions providing the splicing signals for posttranscriptional mRNA processing for expression of target and marker genes in the LEXSY host P10. **SP** designates the signal peptide of *L. mexicana* secreted acid phosphatase LMSAP1 and **H6** the hexa-Histidine stretch. The sequence of the multiple cloning sites is indicated in Figure 2. The DNA Sequence and further description of pLEXSY-hyg2.1 can be downloaded from our website at <https://www.jenabioscience.com/images/103bb272b3/EGE-272.txt>. Please note, that the *KasI* site at the SPCS is not unique and not suited for target gene insertion.




**Figure 5:** Map of **pLEXSY-neo2.1** (Cat. No. EGE-273) expression vector with cloning sites for the target genes replacing the 1 kb stuffer fragment. **5' ssu** and **3' ssu** are regions for homologous recombination into the host chromosome following linearization of the expression plasmid with **SwaI**. **Pr** is a *Leishmania* ribosomal promoter in front of the vector expression cassette. **Utr1** derived from 0.4k-IR of *L. tarentolae aprt*, **utr2** from 1.4k-IR *camCB* and **utr3** from 1.7k-IR are optimized gene-flanking non-translated regions providing the splicing signals for posttranscriptional mRNA processing for expression of target and marker genes in the LEXSY host P10. **SP** designates the signal peptide of *L. mexicana* secreted acid phosphatase LMSAP1 and **H6** the hexa-Histidine stretch. The sequence of the multiple cloning sites is indicated in Figure 2. The DNA Sequence and further description of pLEXSY-neo2.1 can be downloaded from our website at <https://www.jenabioscience.com/images/103bb272b3/EGE-273.txt>. Please note, that the *KasI* site at the SPCS is not unique and not suited for target gene insertion.



**Figure 6:** Map of **pLEXSY-sat2.1** (Cat. No. EGE-274) expression vector with cloning sites for the target genes replacing the 1 kb stuffer fragment. **5'ssu** and **3'ssu** are regions for homologous recombination into the host chromosome following linearization of the expression plasmid with **Swal**. **Pr** is a *Leishmania* ribosomal promoter in front of the vector expression cassette. **Utr1** derived from 0.4k-IR of *L. tarentolae aprt*, **utr2** from 1.4k-IR *camCB* and **utr3** from 1.7k-IR are optimized gene-flanking non-translated regions providing the splicing signals for posttranscriptional mRNA processing for expression of target and marker genes in the LEXSY host P10. **SP** designates the signal peptide of *L. mexicana* secreted acid phosphatase LMSAP1 and **H6** the hexa-Histidine stretch. The sequence of the multiple cloning sites is indicated in Figure 2. The DNA Sequence and further description of pLEXSY-sat2.1 can be downloaded from our website at <https://www.jenabioscience.com/images/103bb272b3/EGE-274.txt>. Please note, that the *KasI* site at the SPCS is not unique and not suited for target gene insertion



**Figure 7:** Map of **pLEXSY-bsd2.1** (Cat. No. EGE-275) expression vector with cloning sites for the target genes replacing the 1 kb stuffer fragment. **5'ssu** and **3'ssu** are regions for homologous recombination into the host chromosome following linearization of the expression plasmid with **SwaI**. **Pr** is a *Leishmania* ribosomal promoter in front of the vector expression cassette. **Utr1** derived from 0.4k-IR of *L. tarentolae aprt*, **utr2** from 1.4k-IR *camCB* and **utr3** from 1.7k-IR are optimized gene-flanking non-translated regions providing the splicing signals for posttranscriptional mRNA processing for expression of target and marker genes in the LEXSY host P10. **SP** designates the signal peptide of *L. mexicana* secreted acid phosphatase LMSAP1 and **H6** the hexa-Histidine stretch. The sequence of the multiple cloning sites is indicated in Figure 2. The DNA Sequence and further description of pLEXSY-bsd2.1 can be downloaded from our website at <https://www.jenabioscience.com/images/103bb272b3/EGE-275.txt>. Please note, that the *KasI* site at the SPCS is not unique and not suited for target gene insertion




**Figure 8:** Map of **pLEXSY-pac2.1** (Cat. No. EGE-276) expression vector with cloning sites for the target genes replacing the 1 kb stuffer fragment. **5' ssu** and **3' ssu** are regions for homologous recombination into the host chromosome following linearization of the expression plasmid with **Swal**. **Pr** is a *Leishmania* ribosomal promoter in front of the vector expression cassette. **Utr1** derived from 0.4k-IR of *L. tarentolae aprt*, **utr2** from 1.4k-IR *camCB* and **utr3** from 1.7k-IR are optimized gene-flanking non-translated regions providing the splicing signals for posttranscriptional mRNA processing for expression of target and marker genes in the LEXSY host P10. **SP** designates the signal peptide of *L. mexicana* secreted acid phosphatase LMSAP1 and **H6** the hexa-Histidine stretch. The sequence of the multiple cloning sites is indicated in Figure 2. The DNA Sequence and further description of pLEXSY-pac2.1 can be downloaded from our website at <https://www.jenabioscience.com/images/103bb272b3/EGE-276.txt>. Please note, that the *KasI* site at the SPCS is not unique and not suited for target gene insertion

## 8.2 Sequences of the primers available for LEXSYcon2.1 kits and illustration of genomic PCR diagnostics

Insert sequencing forward primer <b>P1442</b>	all "AP" expression vectors with 5'utr <i>aprt</i>	5'-CCGACTGCAACAAGGTGTAG-3'	<b>PM-110</b>
Insert sequencing rev. primer <b>A264</b>	all LEXSY expression vectors	5'-CATCTATAGAGAAGTACACGTAAAAG-3'	<b>PM-101</b>
<i>ssu</i> forward primer <b>F3001</b>	5' integration diagnostics of all <i>ssu</i> expression vectors	5'-GATCTGGTTGATTCTGCCAGTAG-3'	<b>PM-105</b>
<i>aprt</i> reverse primer <b>A1715</b>	5' integration diagnostics of all "AP" expression vectors with 5'utr <i>aprt</i>	5'-TATTCGTTGTCAGATGGCGCAC-3'	<b>PM-111</b>
<i>sat</i> forward primer <b>F2999</b>	3' integration diagnostics of all <i>sat</i> expression vectors	5'-CCTAGTATGAAGATTTCCGGTGATC-3'	<b>PM-103</b>
<i>ble</i> forward primer <b>A708</b>	3' integration diagnostics of all <i>ble</i> expression vectors	5'-GGATCCACCGCATGGCCAAGTTGACCAGTG-3'	<b>PM-107</b>
<i>hyg</i> forward primer <b>A3804</b>	3' integration diagnostics of all <i>hyg</i> expression vectors	5'-CCGATGGCTGTGTAGAAGTACTCG-3'	<b>PM-109</b>
<i>neo</i> forward primer <b>H4662</b>	3' integration diagnostics of all <i>neo</i> expression vectors	5'-GGATCCAATATGGGATCGGCCATTG-3'	<b>PM-118</b>
<i>bsd</i> forward primer <b>H4895</b>	3' integration diagnostics of all <i>bsd</i> expression vectors	5'-ATGGCCAAGCCTTTGTCTCAAG-3'	<b>PM-119</b>
<i>pac</i> forward primer <b>H4896</b>	3' integration diagnostics of all <i>pac</i> expression vectors	5'-CACCATGACCGAGTACAAGCC -3'	<b>PM-120</b>
<i>ssu</i> reverse primer <b>F3002</b>	3' integration diagnostics of all <i>ssu</i> integration vectors	5'-CTGCAGGTTACCTACAGCTAC-3'	<b>PM-104</b>

## PCR diagnostics of integration into *ssu* locus of *L. tarentolae*



**Figure 9:** Illustration of genomic PCR diagnostics of *ssu* integration in recombinant LEXSY-2.1 strains. Alternative primers and fragment sizes are shown for the four marker genes. See chapter 5.3 for more details.

### Which primers are included in which kit?

#### LEXSYcon2.1 kit with pLEXSY-sat2.1

- Insert sequencing forward primer **P1442**
- Insert sequencing reverse primer **A264**
- *ssu* forward primer **F3001**
- AP reverse primer **A1715**
- **sat** forward primer **F2999**
- *ssu* reverse primer **F3002**

#### LEXSYcon2.1 kit with pLEXSY-hyg2.1

- Insert sequencing forward primer **P1442**
- Insert sequencing reverse primer **A264**
- *ssu* forward primer **F3001**
- AP reverse primer **A1715**
- **hyg** forward primer **A3804**
- *ssu* reverse primer **F3002**



**LEXSYcon2.1 kit with pLEXSY-ble2.1**

- Insert sequencing forward primer **P1442**
- Insert sequencing reverse primer **A264**
- ssu forward primer **F3001**
- AP reverse primer **A1715**
- **ble** forward primer **A708**
- ssu reverse primer **F3002**

**LEXSYcon2.1 kit with pLEXSY-neo2.1**

- Insert sequencing forward primer **P1442**
- Insert sequencing reverse primer **A264**
- ssu forward primer **F3001**
- AP reverse primer **A1715**
- **neo** forward primer **H4662**
- ssu reverse primer **F3002**

**LEXSYcon2.1 kit with pLEXSY-bsd2.1**

- Insert sequencing forward primer **P1442**
- Insert sequencing reverse primer **A264**
- ssu forward primer **F3001**
- AP reverse primer **A1715**
- **bsd** forward primer **H4895**
- ssu reverse primer **F3002**

**LEXSYcon2.1 kit with pLEXSY-pac2.1**

- Insert sequencing forward primer **P1442**
- Insert sequencing reverse primer **A264**
- ssu forward primer **F3001**
- AP reverse primer **A1715**
- **pac** forward primer **H4896**
- ssu reverse primer **F3002**

### 8.3 Preparation of LEXSY BHI agar plates for clonal selection

For 4 plates prepare 50 ml medium and bring to 37°C

<i>Component</i>	<i>Storage</i>	<i>Amount</i>
<b>2x LEXSY BHI</b> (74 g/L)	room temperature	36 ml
inactivated Fetal Calf Serum (FCS)*	-20°C	10 ml
1M HEPES, pH 7.4	4°C	4 ml
<b>Pen-Strep</b>	-20°C	0.5 ml
<b>Hemin</b> (0,25% in 30% Triethanolamine)	4°C	0.2 ml

Selective antibiotic(s), if necessary

\*Inactivated for 20 min at 56°C (or 1h at 52°C)

- Autoclave or melt (microwave) 50 ml 2% BACTO-Agar (DIFCO) and keep at 55°C
- Pour the medium into the warm agar, mix gently to avoid air bubbles, and distribute 25 ml per plate with a serological pipette (air bubbles may be removed with the pipette)
- Dry the plates after solidifying for 10 min **open** under the laminar flow
- Use the freshly prepared plates immediately, at least on the same day
- Optionally, cover the surface of the plates with nitrocellulose membrane after aspiration.

Refer to chapter 5.1 for clonal selection.

For your convenience, the **LEXSY Plating Kit** (Cat. No. ML-451) is available from Jena Bioscience, containing all components for the preparation of 40 LEXSY BHI agar plates.

## 8.4 Alternative electroporation protocol

### High-Voltage protocol for transfection of LEXSY (after Robinson *et al.* 2003)

- On Friday inoculate *L. tarentolae* pre-culture 1:20 in 10 ml **LEXSY BHI** medium supplemented with Hemin and PenStrep (refer to chapter 3.1) and incubate in tissue culture (TC) flask upright @ 26°C dark until Monday
- On Monday dilute pre-culture 1:10 in 10 ml medium and incubate in TC flask flat @ 26°C o/n
- On Tuesday check cell density of the culture until approx.  $6 \times 10^7$  cells/ml are reached (OD 1.4)\* and ensure by microscopy that the cells are vital and of droplike shape
- Spin cells 3 min, 2000g at room temperature and remove 1/2 volume of supernatant
- Resuspend pellet in remaining medium to get  $10^8$  cells/ml and put on wet ice for 10 min
- Have ready on wet ice in parallel tubes with 1-10 µg of transforming DNA in max. 50 µl of water or Tris buffer pH 8.0 and electroporation cuvettes d=4 mm\*\*
- Add 450 µl pre-chilled cells to the tube with DNA, mix and transfer to the electroporation cuvette on wet ice
- Pulse 2 times at 1500 V, 25 µF with 10 sec between pulses (pulse time ca. 1.2 msec)\*\*\*
- Put cuvette back on ice for exactly 10 min
- Transfer electroporated cells with capillary to 10 ml LEXSY BHI medium in a ventilated TC flask
- Incubate o/n at 26°C as static suspension culture (ca. 20h, OD 0.3-0.4).

\* if the cell density differs from this value, concentrate cells in the next step to get  $10^8$  cells/ml. For transfection, cultures between OD 1.0-1.8 can be used. Do not transfect cells if they are long and thin by microscopy

\*\* use electroporation cuvettes with long electrodes. The entire volume of 0.5 ml must be between the electrodes. Do not use electroporation cuvettes with short electrodes leaving most of the volume outside of the linear electric field


\*\*\* using BioRad GENE PULSER II with PULSE CONTROLLER II and CAPACITANCE EXTENDER PLUS or GENE PULSER Xcell with PC and CE Modules

Low voltage protocol		High voltage protocol
<u>Exptl. decay</u>	<u>Time constant</u>	
d=2mm V=0.4 ml	d=2mm V=0.4 ml	d=4mm V=0.5 ml
set: 450V 450 µF 1 pulse	set: 450V T=3.5ms 1 pulse	set: 1500V 25 µF 2 pulses @ 10" interval
get: 5 - 6 ms ≈20 Ω ≈450 V deliv.	get: ≈3.5 ms ≈20 Ω ≈450 V deliv.	get: ≈1.2 ms ≈50 Ω ≈1500 V delivered

**Figure 10:** Overview of parameters of electroporation protocols with BioRad GENE PULSER. The standard Low voltage electroporation protocol is described in section 4.4


## 8.5 How to grow a *Leishmania* culture

- *L. tarentolae* needs **aerobic** conditions. The strains can be maintained as continuous suspension culture with regular dilutions (refer to chapter 3.2). All culturing is performed at **26°C**. Higher temperatures lower the growth-rates and vitality significantly and *L. tarentolae* will not survive at 37°C in axenic suspension culture
- All growth media should be supplemented with **Hemin** which is essential for *Leishmania*. Hemin is light-sensitive, so *Leishmania* must be cultured in the dark. After completion with Hemin the medium must be stored in the dark at 4°C. For optimal growth and vitality the completed medium should be used within 2 weeks. However, if this shelf life is exceeded, it is possible to re-add Hemin (and PenStrep) and to use this medium for 2 more weeks
- For **maintaining** LEXSY strains for transfection and analysis it is convenient to grow static suspension cultures in 10 ml **LEXSY BHI medium** in ventilated tissue culture (TC) flasks. Don't use agitated cultures for strain maintenance since cells will age much faster. It is not necessary or growth-promoting, to add serum to the BHI medium
- Best results are obtained with inoculations during early stationary phase. Avoid repeated successive dilution of cultures of low cell densities as this may drop growth. However, occasional higher dilutions of stationary cells at e.g. 1:100 will not adversely affect subsequent growth. It is convenient to dilute 10 ml cultures 1:50 on Monday and 1:20 on Friday and to incubate TC flask upright, this way lowering aeration for achieving longer intervals between passages. Don't culture *Leishmania* much longer than for 7 days in the same medium without dilution. For culturing for transfection refer to chapter 4.4, for culturing for protein expression refer to chapter 5.4
- Always control appearance and motility of cells by microscopy. Cells of mid-growth phase cultures are of drop-like shape (Fig. 10), approx. 15x5 µm in size with one flagellum at the flat end, and motile. **These cells are most efficient for transfection and plating on solid media.** Mid-growth phase cultures always contain subpopulations of non- or less motile cells and of cells of different shape. Don't hesitate to transfect, plate or preserve a culture with drop-like cells containing such subpopulations. Cells of older cultures get longer and thinner (needle-like shape) and remain motile. Enhanced motility may result from nutrient deprivation or other limitations and must not necessarily be a sign of mid growth culture stage. Also, bacterial, fungal or other contaminations may be identified by microscopy
- Keep patient, esp. if you are used to working with bacteria. *Leishmania* cells are protozoans with regular doubling times of 7 h in static suspension cultures and 4 - 6 h in agitated cultures. They need their time to grow or to adapt to new conditions
- If you - despite following these instructions - encounter growth problems with the host strain, sediment cells 3 min at 2000g, resuspend pellet carefully in fresh growth medium and continue incubation in ventilated TC flasks. This approach was very helpful in rescuing cultures esp. after transfection
- It is not necessary to centrifuge *Leishmania* cultures at high speed >10.000g. Centrifugation at 2000-3000g is sufficient for sedimentation of cells and makes gentle and quick resuspension of cell pellets easier. However, *Leishmania* cells will survive centrifugation at 20.000g without lysis
- If you culture LEXSY strains in bioreactors be careful with stirring to avoid shear stress. We found it sufficient to aerate the culture in a 10 L fermentation with low speed stirring at 100-300 rpm or even without stirring for obtaining high cell densities up to 10<sup>9</sup> cells/ml. However, *L. tarentolae* was found resistant to shear stress at rotor speeds of more than 500 rpm with turbines consisting of two angular paddles, whereas 700 rpm damaged the cells.



**Figure 11:** Microscopic image of LEXSY cells expressing green fluorescent protein. **A:** exponentially growing culture, suitable for electroporation and cryoconservation **B:** stationary culture, do not use for electroporation or cryopreservation.

## 8.6 Correlation of optical density and cell concentration



**Figure 12:** Calibration of OD readings of spectrophotometer (biowave CO8000) with cell densities determined with a particle counter (Coulter). LEXSY cultures were grown to different stages and OD readings and cell counts were determined from the same sample. The non-linear behaviour at higher OD is due to changes in cell shapes and cell lysis in late stationary cultures.