

## DATA SHEET

**Proteinase K - Solid**

from *Tritirachium album*  
 Endopeptidase K

Cat. Nº.	Amount
<input type="checkbox"/> DPK-103S	100 mg
<input type="checkbox"/> DPK-103M	200 mg
<input type="checkbox"/> DPK-103L	500 mg

**Unit Definition:** One unit is the amount of enzyme which releases at 37 °C in 1 min as many folin-positive amino acids and peptides from haemoglobin as 1 µmol of tyrosine.

**Shipping:**

Shipped on blue ice

**Storage Conditions:**

Store at - 20 °C)

**Shelf life:**

12 months

**Molecular Weight:**

28.9 kDa

**CAS#:**

39450-01-6

**EC number:**

254-457-8

**Purity:**

free of RNases, DNases and Exonucleases

**Form:**

white powder

**Applications:**

Digestion of proteins during DNA and RNA preparation.

**Description:**

Proteinase K is a serine protease that exhibits a very broad cleavage specificity. The protein with a molecular weight of 28.9 kDa cleaves peptide bonds adjacent to the carboxylic group of aliphatic and aromatic amino acids. Proteinase K is not inactivated by metal chelating reagents such as EDTA or detergents such as SDS and is active over a wide range of pH (4 - 12.5).

Proteinase K is a highly active and stable protease with low cutting specificity. The enzyme belongs to the group of subtilisine-related serine proteases and is strongly inhibited by PMSF. In presence of 0.5 - 1 % SDS Proteinase K inactivates DNases and RNases in eucaryotic and microbiological cell cultures. The use of Proteinase K during lysis of the cells allows the isolation of intact highly-molecular nucleic acids.

**Activity:****Specific activity:**

> 30 units/mg

**Dnase activity:**

not detectable

**Rnase activity:**

not detectable