

DATA SHEET



Bacteria DNA Preparation Kit

Spin column based genomic DNA purification from bacteria.

Cat. Nº.	Amount
<input type="checkbox"/> DPK-116S	50 preparations
<input type="checkbox"/> DPK-116L	250 preparations

Shipping:

Shipped at ambient temperature

Storage Conditions:

Store at ambient temperature (except Enzyme Mix and Proteinase K - store at - 20 °C)

Shelf life:

36 months

For in vitro use only!

Kit Contents:

- Cell Resuspension Solution
- Lysis Buffer
- Binding Buffer
- Enzyme Mix (before use, add Cell Resuspension Solution as indicated on the tube and store at -20 °C)
- Proteinase K (before use, add double distilled water as indicated on the tube - store at -20 °C)
- Activation Buffer
- Washing Buffer (before use, add 96-99 % Ethanol as indicated on the tube)
- Elution Buffer
- Spin Columns and 2 ml Collection Tubes

Additional Materials Required:

- 96-99% Ethanol
- Ultra-Pure water
- 1.5 or 2.0 mL microtubes

Applications:

The obtained DNA is suitable for a variety of applications, including PCR amplifications, digestion with restriction enzymes and membrane hybridizations. Before start, add the following components as indicated on the respective bottle/tube:

Not included in the kit:

- double distilled water to Proteinase K
- 96-99 % Ethanol to the Washing Buffer

Included in the kit:

- Cell Resuspension Solution to Enzyme Mix

Description:

The spin column based Bacteria and Plants tissue DNA Preparation Kit is designed for rapid and pure isolation of total DNA from bacteria Gram(+) or Gram(-). The spin column based method completely removes PCR inhibitors such as divalent cations and proteins resulting in a high purity preparation of genomic DNA. There is no use of phenol or chloroform, handling is safe and does not produce any harmful waste.

PREPARATION PROCEDURE:

For S pack (50 preps): Before start, add 500 µl dd-water to the Proteinase K tube, 1,5 ml of Cell resuspension solution to Enzyme Mix tube and 44 ml 96-99 % Ethanol to Washing Buffer bottle.

For L pack (250 preps): Before start, add 500 µl dd-water to each Proteinase K tube, 1,5 ml of Cell resuspension solution to each Enzyme Mix tube and 208 ml 96-99 % Ethanol to Washing Buffer bottle.

Buffer	DPK-116S 50 preps	DPK-116L 250 preps
Lysis Buffer	16 ml	80 ml
Binding Buffer	16 ml	80 ml
Enzyme Mix	2 microtubes	10 microtubes
Proteinase K (10 mg/mL)	5 mg	5 x 5 mg
Cell Resuspension buffer	3,2 ml	16 ml
Activation Buffer	6 ml	30 ml
Washing Buffer	11 ml Add 44ml EtOH	52 ml Add 208 ml EtOH
Elution Buffer	5 mL	25 mL

DATA SHEET



Bacteria DNA Preparation Kit

Spin column based genomic DNA purification from bacteria.

1. PREPARATION FROM BACTERIA

It is essential to use the correct amount of starting material in order to obtain optimal DNA yield and purify. A maximum amount of 10^9 bacteria cells can generally be processed. Overnight cultured bacteria cells can be processed. Cell pellets can be stored at -70°C for several months.

A. Cell Lysis

- Harvest maximum 1 ml of cultured bacteria cells by centrifugation at 10,000 g for 1 min
- Discard the supernatant and vigorously vortex.

NOTE: After discarding the supernatant, the remaining 10-50 μl of media broth doesn't affect the next purification step. The vigorous vortexing will resuspend the pelleted cells in the remaining media and help the cell lysis

- Add **50 μl of Enzyme Mix** and vortex vigorously for 30-60 sec.
- Incubate tube at 37°C for 10 min.
- Add **300 μl of Lysis buffer**
- Add **8 μl of Proteinase K** and mix by pipetting.
- Incubate at 60°C for 10 min and cool down on ice for 5 min.

B. DNA Binding

- Add **300 μl Binding Buffer** and vortex briefly or mix by inverting
- Place the tube on ice for 5 min
- Centrifuge for 5 min at 10,000 g (repeat if the supernatant is not clear)

C. Column Activation

- Place a spin column into a 2 ml collection tube
- Add **100 μl Activation Buffer** into the Spin Column
- Centrifuge at 10,000 g for 30 sec and immediately proceed to next step
- Discard the flow-through

D. Column Loading:

- Transfer 600 μl of the supernatant from **section B. DNA Binding** into a clean 1.5 ml tube.
- Add 200 μl of absolute ethanol and vortex vigorously.
- Pipet **400 μl of the mixture** directly into the Spin Column (to avoid overloading of the column).
- Centrifuge for 1 min at 10,000 g.
- Discard the flow-through.
- Pipet the remaining **400 μl of the mixture** into the Spin Column.
- Centrifuge for 1 min at 10,000 g.
- Discard the flow-through.

E. Washing:

- Add **500 μl Washing Buffer** into spin column
- Centrifuge for 30 sec at 10,000 g
- Discard the flow-through

Optional: Repeat Washing Step if highest purity is required

F. Remove residual Washing Buffer

- Centrifuge at 10,000 g for 2 min to remove residual Washing Buffer
- Discard the 2 ml Collection Tube
- Place the Spin Column into a new 1.5 ml microtube

G. Elution of DNA:

- Add **40-50 μl Elution Buffer** into the center of the spin column
- Incubate at room temperature for 1 min
- Centrifuge at 10,000 g for 2 min. Discard the spin column
- Store DNA at 4°C or -20°C

DATA SHEET



Bacteria DNA Preparation Kit

Spin column based genomic DNA purification from bacteria.

Troubleshooting

Problem	Cause	Suggestions
Low yield	Insufficient lysis	Reduce the amount of starting material. Prolong the lysis step. Shake the sample during incubation
	No Ethanol added to Washing Buffer	Make sure to add EtOH to the Washing Buffer before start
	DNA is sheared/degraded	Avoid excessive pipetting after lysis. Do not vortex after lysis
	DNA is not completely eluted	Elute in two steps. Increase the volume of Elution Buffer. Prolong the incubation of the Elution Buffer on the Spin Column to 5 min. Make sure the Elution Buffer is delivered onto the matrix and not on the wall of the Spin Column.
RNA contamination	No / not enough Enzyme Mix added	Make sure to add Enzyme Mix in the Cell Lysis step
Inhibition of downstream enzymatic reactions	Ethanol carryover	Make sure to remove residual Washing Buffer before elution. Remove the Spin Column carefully from the Collection Tube so that it does not come in contact with the flow-through. If in doubt, centrifuge again for 1 min.