

DATA SHEET

**Total RNA Purification Kit**

Isolation of total RNA by silica-gel membrane adsorption

Cat. Nº.	Amount
<input type="checkbox"/> DPK-108XS	10 preparations
<input checked="" type="checkbox"/> DPK-108S	50 preparations
<input type="checkbox"/> DPK-108L	250 preparations

Shipping:

Shipped at ambient temperature

Storage Conditions:

Store at ambient temperature

Shelf life:

12 months

For *in vitro* use only!**Kit Contents:**

- Lysis Buffer (before use, add 2-Mercaptoethanol as indicated on the bottle) - stable for 1 month at room temperature
- Activation Buffer
- Primary Washing Buffer (before use, add 96-99 % Ethanol as indicated on the bottle)
- Secondary Washing Buffer (before use, add 96-99 % Ethanol as indicated on the bottle)
- Elution Buffer
- Spin Columns
- 2 ml Collection Tubes

Additional Materials Required:

- 2-Mercaptoethanol (2-ME)
- Optional: Chloroform
- 96-99 % Ethanol
- 2-Propanol (Isopropanol)
- 1.5 ml microtubes

Description:

Total RNA Purification Kit is designed for rapid, high purity and high yield isolation of total RNA from small amounts of various samples including blood, animal and plant tissue, bacteria and viruses. The spin column based method allows complete removal of inhibitors such as divalent cations and proteins. Due to elimination of phenol, handling of the kit is safe and no harmful waste is produced. The purified total RNA can be used in a number of downstream applications. The kit allows the purification of up to 100 µg RNA per preparation.

Preparation Procedure:

Before start, add the following components (not included in the kit) as indicated on the respective bottle:

- 2-Mercaptoethanol to the Lysis Buffer (10 µl 2-Mercaptoethanol per 1 ml Lysis Buffer)
- 96-99 % Ethanol to Primary Washing Buffer and Secondary Washing Buffer

Buffer	DPK-108XS 10 preps	DPK-108S 50 preps	DPK-108L 250 preps
Lysis Buffer	5.2 ml (add 52 µl 2-ME)	26 ml (add 260 µl 2-ME)	130 ml (add 1.3 ml 2-ME)
Activation Buffer	1.2 ml	6 ml	30 ml
Primary Washing Buffer	add 1.6 ml Ethanol (final volume 8 ml)	add 8 ml Ethanol (final volume 40 ml)	add 40 ml Ethanol (final volume 200 ml)
Secondary Washing Buffer	add 6.4 ml Ethanol (final volume 8 ml)	add 32 ml Ethanol (final volume 40 ml)	add 160 ml Ethanol (final volume 200 ml)
Elution Buffer	1 ml	5 ml	25 ml

1. Sample Preparation and Cell Lysis:Fresh Tissue Sample - Animals or Plants

- Collect 20-50 mg fresh tissue sample in a micro-centrifuge tube.
- Add 300 µl of Lysis Buffer (2-ME added) and homogenize the material using an appropriate apparatus (hand-operated pellet pestle or motor-driven grinder).
- Add additional 200 µl of Lysis Buffer (2-ME added) to the homogenized sample and vortex 15-30 s (Note: Sample volume should not exceed 10 % of the Lysis Buffer volume).
- Centrifuge at 10,000 g for 10 min.

DATA SHEET



Total RNA Purification Kit

Isolation of total RNA by silica-gel membrane adsorption

Optional step in case that debris still remains in the supernatant:

- Add 500 µl chloroform (not included in the kit) and vortex for 15-30 s.
- Centrifuge at 10,000 g for 10 min.
- Transfer the supernatant (if you added chloroform: the upper aqueous phase) into a microcentrifuge tube.

Blood

- Transfer 100 µl of non-coagulating blood to a microcentrifuge tube.
- Add 500 µl of Lysis Buffer (2-ME added) and vortex for 10 s.

Cells from Nasal or Throat Swabs

- Add 500 µl of Lysis Buffer (2-ME added) to a microcentrifuge tube.
- Brush a sterile, single-use cotton swab or Buccal Swab Brush inside the nose or mouth of the subject.
- Cut the cotton tip where the nasal or throat cells were collected and place it into the microcentrifuge tube containing the Lysis Buffer (2-ME added).
- Close the tube, vortex and incubate at room temperature for 5 min.

Cells Grown in Monolayer

- Put off culture media.
- Add 500 µl of Lysis Buffer (2-ME added) per $1-5 \times 10^6$ cells.
- Lyse cells and homogenize the sample by pipetting up and down several times.

Cells Grown in Suspension

- Pellet $1-5 \times 10^6$ animal, plant or yeast cells, or 1×10^7 bacterial cells. (Occasionally, enzymatic lysis or mechanical disruption is required for cell-wall disruption of some yeast and bacterial cells.)
- Discard the supernatant and add 500 µl of Lysis Buffer (2-ME added).
- Lyse the sample by repetitive pipetting or vortexing for 10 s.

2. Column Activation [optional]

- Place a spin column into a 2 ml collection tube.
- Add 100 µl Activation Buffer into the Spin Column.
- Centrifuge at 10,000 g for 30 s.
- Discard the flow-through.

3. Column Loading:

- Add 300 µl (or 0.6 x volume of the cell lysate) Isopropanol to the prepared lysate and vortex.
- Transfer the mixture directly into the spin column.
- Centrifuge at 10,000 g for 30 s.
- Discard the flow-through.

4. Primary Column Washing:

- Apply 700 µl of Primary Washing Buffer (ethanol added) to the Spin Column.
- Centrifuge at 10,000 g for 30 s.
- Discard the flow-through.

5. Secondary Column Washing:

- Apply 700 µl of Secondary Washing Buffer (ethanol added) to the Spin Column.
- Centrifuge at 10,000 g for 30 s.
- Discard the flow-through.
- Centrifuge again at 10,000 g for 2 min to remove residual ethanol.

6. Elution of RNA:

- Place the Spin Column into a DNase/RNase-free microcentrifuge tube.
- Add 40-50 µl Elution Buffer to the center of the column membrane.
- Incubate at room temperature for 1 min.
- Centrifuge at 10,000 g for 1 min to elute the RNA
- Store RNA at -20 or -80 °C.

Elimination of remaining DNA:

Remaining genomic DNA may be particularly a problem in subsequent RT-PCR or quantification of low-copy transcripts. For complete removal of gDNA from RNA preparations **gDNA Removal Kit (Cat.-No. PP-219 from Jena Bioscience)** is recommended. The kit is based on a heat labile dsDNase which is irreversibly inactivated at moderate temperatures.