

DATA SHEET



Blood-Animal-Plant DNA Preparation Kit

Spin column based genomic DNA purification from blood, animal and plant cells

Cat. Nº.	Amount
<input type="checkbox"/> DPK-112S	50 preparations
<input type="checkbox"/> DPK-112L	250 preparations

Shipping:

Shipped at ambient temperature

Storage Conditions:

Store at ambient temperature (except RNase and Proteinase K - store at - 20 °C)

Shelf life:

12 months

For *in vitro* use only!

Kit Contents:

- Blood Lysis Buffer
- Lysis Buffer
- Binding Buffer
- Rnase A (before use, add pure water as indicated on the bottle) - store at -20°C
- Proteinase K (before use, add pure water as indicated on the bottle) - store at -20°C
- Activation Buffer
- Washing Buffer (before use, add 96-99% Ethanol as indicated on the bottle)
- Elution Buffer
- Spin columns and 2 mL Collection Tubes

Additional Materials Required:

- 96-99% Ethanol
- Pure water
- 1.5 or 2.0 mL microtubes

Applications:

The obtained DNA is suitable for a variety of applications, including real-time PCR, southern blot analysis, genotyping and discovery or validation of SNP/SSR markers.

Before start, add the following components (not included in the kit) as indicated on the respective bottle/tube:

- Pure water to RNase A and Proteinase K
- 96-99% Ethanol to the Washing Buffer

Description:

The spin column based Blood-Animal-Plant DNA Preparation Kit is designed for rapid and high purity isolation of genomic DNA from whole blood, animal cells and plant tissue. The spin column based method completely removes PCR inhibitors such as divalent cations and proteins resulting in a high purity preparation of genomic DNA. There is no use of phenol or chloroform, handling is safe and does not produce any harmful waste. Column based genomic DNA purification kits yield up to 30 µg DNA.

Preparation Procedure:

For XS pack (10 preps): Before start, add 500 µL pure water to the Proteinase K tube, 150 µL pure water to the RNase A tube and 8 mL 96-99% Ethanol (not included in the kit) to the Washing Buffer bottle.

For S pack (50 preps): Before start, add 500 µL pure water to the Proteinase K tube, 150 µL pure water to the RNase A tube and 40 L 96-99% Ethanol (not included in the kit) to the Washing Buffer bottle.

For L pack (250 preps): Before start, add 500 µL pure water to each Proteinase K tube, 150 µL dd-water to each RNase A tube and 200 mL 96-99% Ethanol (not included in the kit) to each Washing Buffer bottle.

Buffer	DPK-112XS 10 preps	DPK-112S 50 preps	DPK-112L 250 preps
Blood Lysis Buffer	11 mL	55 mL	260 mL
Lysis Buffer	3.2 mL	16 mL	80 mL
Binding Buffer	3.2 mL	16 mL	80 mL
Rnase A (50 mg/mL)	7.5 mg	7.5 mg	5 x 7.5 mg
Proteinase K (10 mg/mL)	5 mg	5 mg	5 x 5 mg
Activation Buffer	1.1 mL	5.5 mL	27 mL
Washing Buffer	add 8 mL Ethanol (final volume 10 mL)	add 40 mL Ethanol (final volume 50 mL)	add 200 mL Ethanol to each bottle (final volume 250 mL)
Elution Buffer	0.55 mL	2.7 mL	13.5 mL

1. DNA Preparation from blood

A. Sample Preparation

- Mix 200 µL (1 volume) of whole human blood with 1 mL (5 volumes) of Blood Lysis Buffer in an appropriately sized tube
- Incubate for 10 min on ice
- Mix by vortexing briefly 2-3 times during incubation
- Centrifuge at 10,000 g for 10 min to pellet the white blood cells
- Remove or discard supernatant completely

DATA SHEET



Blood-Animal-Plant DNA Preparation Kit

Spin column based genomic DNA purification from blood, animal and plant cells

B. Cell Lysis

- Add 300 μ L Lysis Buffer and 2 μ L RNase A to the cell pellet
- Vortex vigorously for 30-60 s
- Add 8 μ L Proteinase K and mix by pipetting
- Incubate at 60 °C for 10 min and cool down for 5 min
- Add 300 μ L Binding Buffer and vortex briefly
- Place tube on ice for 5 min
- Centrifuge for 5 min at 10,000 g

C. Column Activation [optional]

- Place a spin column into a 2 mL collection tube
- Add 100 μ L Activation Buffer into the spin column
- Centrifuge at 10,000 g for 30 s and immediately proceed to the next step
- Discard flow-through

D. Column Loading:

- Pipet supernatant directly into the spin column
- Centrifuge for 1 min at 10,000 g
- Discard flow-through

E. Primary Washing:

- Add 500 μ L Washing Buffer into the spin column
- Centrifuge for 30 s at 10,000 g
- Discard flow-through

F. Secondary Washing:

- Add 500 μ L Washing Buffer into the spin column
- Centrifuge for 30 s at 10,000 g
- Discard flow-through
- Centrifuge again at 10,000 g for 1 min to remove residual Washing Buffer
- Discard the 2 mL wash tube and place column into the elution tube (new tube)

G. Elution of DNA:

- Add 40-50 μ L Elution Buffer into the center of the column
- Incubate at room temperature for 1 min
- Centrifuge at 10,000 g for 2 min
- Store the eluted DNA at 4 °C or -20 °C

2. DNA Preparation from animal cells

It is essential to use the correct amount of starting material in order to obtain optimal DNA yield and purity. A maximum amount of 80 mg animal material can generally be processed. We recommend to start with no more than 80 mg fresh or frozen animal tissue. Tissues can be stored at -70 °C for several months. Quick-freeze tissues in liquid nitrogen and grind thoroughly with a mortar and pestle. Transfer tissue powder into a 1.5 mL microcentrifuge tube cooled with liquid nitrogen.

A. Cell Lysis:

- Add 300 μ L Lysis Buffer and 2 μ L RNase A to the tissue powder
- Vortex vigorously for 30-60 s
- Add 8 μ L Proteinase K and mix by pipetting
- Incubate at 60 °C for 20 min and cool down for 5 min
- Add 300 μ L Binding Buffer and vortex briefly
- Place the tube one ice for 5 min
- Centrifuge for 5 min at 10,000 g

B. Column Activation [optional]:

- Place a spin column into a 2 ml collection tube
- Add 100 μ L Activation Buffer into the Spin Column
- Centrifuge at 10,000 g for 30 s and immediately proceed to the next step
- Discard flow-through

C. Column Loading:

- Pipet supernatant directly into the spin column
- Centrifuge for 1 min at 10,000 g
- Discard flow-through

D. Primary Washing:

- Add 500 μ L Washing Buffer into spin column
- Centrifuge for 30 s at 10,000 g
- Discard flow-through

E. Secondary Washing:

- Add 500 μ L Washing Buffer into the spin column
- Centrifuge for 30 s at 10,000 g
- Discard flow-through
- Centrifuge again at 10,000 g for 1 min to remove residual Washing Buffer
- Discard the 2 mL wash tube and place column into the elution tube (new tube)

F. Elution of DNA

- Add 40-50 μ L Elution Buffer into the center of the column
- Incubate at room temperature for 1 min
- Centrifuge at 10,000 g for 2 min
- Store the eluted DNA at 4 °C or -20 °C

3. DNA Preparation from plant cells

It is essential to use the correct amount of starting material in order to obtain optimal DNA yield and purity. A maximum amount of 80 mg plant material can generally be processed. We recommend to start with no more than 80 mg fresh or frozen plant tissue. Tissues can be stored at -70 °C for several months. Quick-freeze tissues in liquid nitrogen and grind thoroughly with a mortar and pestle. Transfer tissue powder into a 1.5 mL microcentrifuge tube cooled with liquid nitrogen.

DATA SHEET



Blood-Animal-Plant DNA Preparation Kit

Spin column based genomic DNA purification from blood, animal and plant cells

A. Cell Lysis:

- Add 300 μ L Lysis Buffer and 2 μ L RNase A to the tissue powder
- Vortex vigorously for 30-60 s
- Add 8 μ L Proteinase K and mix by pipetting
- Incubate at 60 °C for 20 min and cool down for 5 min
- Add 300 μ L Binding Buffer and vortex briefly
- Place the tube one ice for 5 min
- Centrifuge for 5 min at 10,000 g

B. Column Activation [optional]:

- Place a spin column into a 2 mL collection tube
- Add 100 μ L Activation Buffer into the Spin Column
- Centrifuge at 10,000 g for 30 s and immediately proceed to the next step
- Discard flow-through

C. Column Loading:

- Pipet supernatant directly into the spin column
- Centrifuge for 1 min at 10,000 g
- Discard flow-through

D. Primary Washing:

- Add 500 μ L Washing Buffer into spin column
- Centrifuge for 30 s at 10,000 g
- Discard flow-through

E. Secondary Washing:

- Add 500 μ L Washing Buffer into the spin column
- Centrifuge for 30 s at 10,000 g
- Discard flow-through
- Centrifuge again at 10,000 g for 1 min to remove residual Washing Buffer
- Discard the 2 mL wash tube and place column into the elution tube (new tube)

F. Elution of DNA

- Add 40-50 μ L Elution Buffer into the center of the column
- Incubate at room temperature for 1 min
- Centrifuge at 10,000 g for 2 min
- Store the eluted DNA at 4 °C or -20 °C