

## DATA SHEET



## Bst 2.0 Polymerase

Bst polymerase optimized for Loop-Mediated Isothermal DNA Amplification (LAMP)

Cat. Nº.	Amount
<input type="checkbox"/> POL-133XS	500 units
<input type="checkbox"/> POL-133S	2.000 units
<input type="checkbox"/> POL-133L	10.000 units

**Concentration:**

8 units/ $\mu$ L

**Shipping:**

Shipped on blue ice

**Storage Conditions:**

Store at -20 °C

**For *in vitro* use only!**

**Additional Storage Conditions:**

Avoid freeze/thaw cycles

**Shelf Life:**

12 months

**Description:**

Bst 2.0 Polymerase is a genetically enhanced Bst polymerase of the next generation. The polymerase is the ideal choice for robust amplification of DNA at constant temperature (60 to 65 °C). The enzyme shows high strand displacement activity and generates an amplification factor of up to  $10^9$  which is comparable to approx. 30 cycles in a PCR assay. The polymerase allows detection of a target gene within 30-60 minutes.

**Kit contents:**
**Bst 2.0 Polymerase (blue cap)**

8 units/ $\mu$ L Bst DNA Polymerase in 10 mM Tris-HCl, 50 mM KCl, 50 % (v/v) Glycerol, pH 7.5 (25 °C) and stabilizers.

**Isothermal Buffer (red cap) - 10x conc.**

200 mM Tris-HCl pH 8.8, KCl,  $(\text{NH}_4)_2\text{SO}_4$ , 20 mM  $\text{MgSO}_4$  and detergents.

 **$\text{MgSO}_4$  Stock Solution (yellow cap)**

25 mM  $\text{MgSO}_4$ .

**Detection**

Although some methods have been developed to visualize DNA amplification by basic equipment or even the naked eye (increase of turbidity, color change of added dyes, hybridization to gold-bound ss-DNA) in general real-time detection of the DNA amplification by a fluorescent DNA-intercalator dye is recommended. Addition of EvaGreen Fluorescent DNA Stain (#PCK-122, not provided) to the assay allows a sensitive measurement of the increasing amount of DNA without influence on the reaction.

**Assay design**

Isothermal amplification is an extremely sensitive detection method and care should be taken to avoid contamination of set-up areas and equipment with DNA of previous reactions. A common problem is amplification in no-template controls due to:

1. carry-over contamination or
2. amplification of unspecifically annealed primers or primer dimer formations.

As sensitivity and non-template amplification of in-silico designed primers may vary, the evaluation of 2-4 real primer sets before choosing a final set is recommended.

## DATA SHEET



## Bst 2.0 Polymerase

Bst polymerase optimized for Loop-Mediated Isothermal DNA Amplification (LAMP)

### Assay set-up

Depending on the detection method and machine a reaction volume of 20-50  $\mu\text{l}$  is recommended for most applications. Pipet with sterile filter tips and perform the set-up in an area separate from DNA preparation or analysis. No-template controls should be included in all amplifications. First, prepare a 10x conc. primer pre-mix. A 10x primer mix should contain: 16  $\mu\text{l}$  FIP, 16  $\mu\text{M}$  BIP, 2  $\mu\text{M}$  F3, 2  $\mu\text{M}$  B3, 4  $\mu\text{M}$  LoopF, 4  $\mu\text{M}$  Loop B in TE buffer or water. Second, set-up the isothermal amplification assay:

Components	25 $\mu\text{l}$ reaction	[ final ]
Isothermal Buffer (10x)	2,5 $\mu\text{l}$	1x
MgSO <sub>4</sub> (25 mM)	4 - 6 $\mu\text{l}$	6 - 8 mM
dNTP (Mix 10 mM)	3,5 $\mu\text{l}$	1,4 mM
Primer Mix (10x)	1,25 - 2,5 $\mu\text{l}$	0,5 - 1x
Bst 2.0 (8 U/ $\mu\text{l}$ ) <sup>*</sup>	1 $\mu\text{l}$	0,32 units/ $\mu\text{l}$
EvaGreen DNA Stain 50x	0,5 $\mu\text{l}$	1 x
Template DNA	x $\mu\text{l}$	1 pg - 10 ng
PCR Grade water	fill up to 25 $\mu\text{l}$	--

\* May vary depending on the application.

- Use a specific detection instrument for isothermal amplification or a real-time PCR cycler to run the assays
- Set the instrument to a constant incubation temperature between 60 to 65°C (depending on the primer annealing temperature)
- Measure the fluorescence intensity at an interval of 1 min for up to 60 min

### Optimization of MgSO<sub>4</sub> concentration:

Isothermal buffer provides a final Mg<sup>2+</sup> concentration of 2 mM. However, most LAMP reactions requires a final Mg<sup>2+</sup> concentration between 6 - 8 mM. To optimize Mg<sup>2+</sup>, add 25 mM of MgSO<sub>4</sub> stock solution as shown in the table.

final MgSO <sub>4</sub> conc.	25 $\mu\text{L}$ reaction
6 mM	4 $\mu\text{l}$
7 mM	5 $\mu\text{l}$
8 mM	6 $\mu\text{l}$

### Trouble shouting

If amplification in no-template controls occurs the following points should be reviewed.

#### Cross contamination from environments

- Clean equipment and areas with "DNA Away" solution
- Replace reagent stocks and pre-mixes with new components
- Stop reactions at an earlier point of time before non-template amplification occur

#### Carry-over contamination from previous reaction products

- Avoid opening reaction vessels after amplification
- Use separate preparation area and equipment if post-reaction processing is necessary
- Add UNG Heat Labile (#ENZ-126, not provided) in your reaction.

#### Non-template amplification from primers

- Increase incubation temperature stepwise by 1-2 °C
- Design a new set of primers for the target sequence