

## DATA SHEET



## FastPol HF DNA Polymerase

High Fidelity DNA polymerase with enhanced processivity.

*Pyrococcus furiosus*, recombinant, *E. coli*

Cat. Nº.	Amount
<input type="checkbox"/> POL-132XS	100 units
<input type="checkbox"/> POL-132S	250 units
<input checked="" type="checkbox"/> POL-132M	500 units
<input type="checkbox"/> POL-132L	2 x 500 units
<input type="checkbox"/> POL-132XL	4 x 500 units

### Unit Definition:

One unit is defined as the amount of the enzyme required to catalyze the incorporation of 10 nmoles of dNTPs into an acid-insoluble form in 30 minutes at 70 °C.

### Concentration:

2 units/ $\mu$ L

### Shipping:

Shipped on blue ice

### Storage Conditions:

Store at -20 °C

**For *in vitro* use only!**

### Additional Storage Conditions:

Avoid freeze/thaw cycles

### Shelf Life:

24 months

### Kit contents:

#### FastPol HF DNA Polymerase (blue cap)

2 units/ $\mu$ L FastPol DNA Polymerase in Tris-HCl, KCl, EDTA, DTT, 50% (v/v) Glycerol, pH 8.0 (25°C) and stabilizers.

#### FastPol HF Reaction Buffer complete (red cap) - 5x conc.

Optimized buffer for FastPol Polymerase.

### Description:

FastPol HF is a designed highly thermostable DNA Polymerase that combines maximum fidelity and processivity in one. With an error rate 6x lower than that of *Pyrococcus furiosus* (Pfu) and extension rate up to 50x greater than Taq, **FastPol** HF DNA Polymerase generates improved product yields with high speed without compromising accuracy. The combined enhanced fidelity and processivity with high yields using minimum amount of enzyme make from FastPol an ideal choice for routine PCR, cloning and also long and difficult amplifications. FastPol HF is supplied with an optimized 5x buffer system containing  $Mg^{2+}$  and suitable for most applications.

\* encoded by the same gene of Phusion®

### Properties:

5'  $\rightarrow$  3' Exonuclease: No

3'  $\rightarrow$  5' Exonuclease: Yes

Product overhang: blunt

Recommended extension time: 15-30 sec/kb

GC-Rich samples: Yes

Long range PCR: up to 20 kbp

### PCR Reaction Setup

The PCR procedure below shows appropriate volumes for a single 50- $\mu$ L reaction. For multiple reactions, prepare a master mix of components common to all and then dispense appropriate volumes into each PCR reaction tube prior to adding template DNA and primers.

Thaw, mix, and briefly centrifuge each component before use.

Add the following components to a microcentrifuge tube:

#### 1. Prepare PCR master mix

Note: Consider the volumes for all components listed next steps to determine the correct amount of water required to reach your final reaction volume.

Components	50 $\mu$ L rxn	[ final ]
Water, grade PCR	To 50 $\mu$ L	
5x Reaction Buffer	10 $\mu$ L	1X
dNTP (Mix 10 mM)	1 $\mu$ L	200 $\mu$ M each
FastPol HF DNA Polymerase (2U/ $\mu$ L)	0,5 $\mu$ L	1 U/reaction

\*Do not exceed 1U/50  $\mu$ L reaction of FastPol

Mix and briefly centrifuge the components.

#### 2. Add template DNA and primers

Components	50 $\mu$ L rxn	[ final ]
Forward primer (10 $\mu$ M)	1,5 - 5 $\mu$ L	0,3 - 1 $\mu$ M
Reverse primer (10 $\mu$ M)	1,5 - 5 $\mu$ L	0,3 - 1 $\mu$ M
DNA template recommended		10 ng**

\* The minimum recommended primer concentration is 0,3  $\mu$ M . For maximum product yield a final concentration of 1  $\mu$ M is used

\*\*Genomic DNA: 10 ng - 1  $\mu$ g; plasmidial or viral DNA: 5 pg - 10 ng

## DATA SHEET

**FastPol HF DNA Polymerase**

High Fidelity DNA polymerase with enhanced processivity.

*Pyrococcus furiosus*, recombinant, *E. coli*

**3. Cycling instruction**

Recommended cycling conditions:

Step	Temp.	Time
Initial denaturation	98 °C	2 min
25 - 35 cycles	Denaturation	98 °C 20 sec
	Annealing <sup>1</sup>	49-68 °C 15 sec
	Elongation	68 °C 15-30 sec/kb
Final extension	68 °C	1 - 10 min
Hold	4 - 8 °C	hold

<sup>1</sup> The annealing temperature depends on the melting temperature of the primers used.

For optimal specificity and amplification an individual optimization of the recommended parameters may be necessary for each new template DNA and/or primer pair.