

N¹-Methylpseudo-UTP NU-890-THR

Product Name

Long Name
Catalog Number
LOT Number
Formula
Molecular Weight
Exact Mass
Spectroscopic Properties
Storage
Shelf Life

N¹-Methylpseudo-UTP

N¹-Methylpseudouridine-5'-triphosphate, Sodium Salt
NU-890-THR
XXX
C₁₀H₁₇N₂O₁₅P₃ (free acid)
498.17 g x mol⁻¹ (free acid)
497.98 g x mol⁻¹ (free acid)
λ_{max} = 271 nm; ε = 8.9 L x mmol⁻¹ x cm⁻¹ (Tris-HCl, pH 7.5)
store at -20 °C
12 months after date of delivery

Test

Form
Color
Concentration
pH
Purity
Identification
Functionality (T7 IVT, 1400 nt
transcript, fluorescence microplate
assay)
DNAses (gel electrophoresis)
RNAses (gel electrophoresis)
Protease (UV-Vis)
Endotoxin (LAL test)

Specification

solution in water
colorless to slightly yellow
100 mM – 105 mM
7.5 ± 0.25
≥ 99 % (RP-HPLC)
MS negative and/or positive mode
≥ 100 µg RNA (100% substitution)
not detectable (LOD: 0.002 U DNase I equiv.)
not detectable (LOD: 1 pg RNase A equiv.)
not detectable (LOD: 5 ng Proteinase K equiv.)
≤ 1.0 EU/ml

Result

conforms
conforms
conforms
conforms
XY %
conforms
conforms
conforms
conforms
conforms