

## DATA SHEET

**Tris - EDTA Buffer pH 7.6 - 1X**

Molecular biology grade

Cat. Nº.	Amount
<input type="checkbox"/> BUF-108S	100 mL
<input type="checkbox"/> BUF-108L	1 L

**Shipping:**

Shipped at ambient temperature

**Storage Conditions:**

Store at ambient temperature

**Concentration:**

1 X

**Shelf life:**

12 months

**Description:**

The 1 x concentration of Tris-EDTA Buffer contains:

10 mM Tris

1 mM EDTA

Tris-EDTA Buffer is commonly used as storage buffer for DNA and RNA samples. The Buffer protects DNA and RNA from degradation.

**Composition:**

10 mM Tris-HCl pH 7.6 and 1 mM EDTA.