

**Ammonium Acetate Solution 5 M**

Cat. No.	Amount
BU-101	100 ml

For general laboratory use.

Shipping: shipped at ambient temperature

Storage Conditions: store at 4 °C

Shelf Life: 12 months

Molecular Formula: CH₃CO₂NH₄

Molecular Weight: 77.08 g/mol

CAS#: 631-61-8

EC number: 211-162-9

Form: liquid

Concentration: 5 M

Description:

Ammonium Acetate is recommended for applications like DNA purification, precipitation or protein crystallization. It is particularly suitable for removal of nucleotides.

Quality:

- Molecular biology grade
- sterile filtered, filter pore size 0.22 µm
- DNase / RNase free