

**Sodium Chloride Solution 5 M**

Cat. No.	Amount
BU-114-5M	100 ml

For general laboratory use.

Shipping: shipped at ambient temperature

Storage Conditions: store at 4 °C

Shelf Life: 12 months

Molecular Formula: NaCl

Molecular Weight: 58.44 g/mol

CAS#: 7647-14-5

Form: liquid

Concentration: 5 M

Description:

Sodium Chloride is used for common cell culture procedures. Applications of NaCl include the removal of small nucleic acid fragments from plasmid DNA preparations and the precipitation of DNA from SDS-containing samples.

Quality:

- Molecular biology grade
- sterile filtered, filter pore size 0.22 µm
- DNase / RNase free