


**AF555-Picolyl-Azide**

Abs/Em = 555/572 nm
structural analog to Alexa Fluor® 555

| Cat. No. | Amount |
|-----------|----------|
| CLK-091-1 | 1 mg |
| CLK-091-5 | 5 x 1 mg |



excitation and emission spectrum of AF555

For general laboratory use.

Shipping: shipped at ambient temperature

Storage Conditions: store at -20 °C

Shelf Life: 12 months after date of delivery

Molecular Weight: 1049.21 g/mol (protonated)

Purity: ≥ 90 % (HPLC)

Form: solid

Solubility: DMSO, DMF

Spectroscopic Properties: λ_{exc} 555 nm, λ_{em} 572 nm,
 ϵ 155.0 L mmol⁻¹ cm⁻¹