

**LEXSY in vitro Translation Cell Extract**

for cell-free protein synthesis

Cat. No.	Amount
EGE-260	60 reactions (2 x 500 µl)

For general laboratory use. Not intended for human or animal diagnostic or therapeutic uses.

Shipping: shipped on dry ice

Storage Conditions: store at -80 °C

Shelf Life: 12 months

Description:

The LEXSY cell extract for *in vitro* translation contains functional ribosomes and all essential components of the eukaryotic translation and folding machinery. For target mRNA generation the heterologous T7 RNA polymerase was added to the extracts.

Please note:

EGE-260 is part of the LEXSY *in vitro* Translation Kit available as plasmid based (EGE-2002-15) or PCR based (EGE-2010-15) cell-free protein synthesis kit.