




pLEXY-ble2.1

Integrative constitutive expression vector
Antibiotic selection of transfectants with LEXSY Bleo

Cat. No.	Amount
EGE-271	5 µg

Description:

General expression vector pLEXY-ble2.1 with the *bleomycin resistance* gene from *S. hindustans* allowing selection of recombinant LEXSY strains with the antibiotic LEXSY Bleo (Cat.-No. AB-103).



For general laboratory use. Not intended for human or animal diagnostic or therapeutic uses.

Shipping: shipped on gel packs

Storage Conditions: Store at -20 °C

Shelf Life: 12 months