




## pLEXY-hyg2.1

integrative constitutive expression vector  
antibiotic selection of transfectants with LEXSY Hygro

Cat. No.	Amount
EGE-272	5 µg

### Description:

General expression vector pLEXY-hyg2.1 with the *hyg* marker gene (encoding Hygromycin phosphotransferase, HPT) allowing selection of recombinant LEXSY strains with the antibiotic LEXSY Hygro (Cat.-No. AB-104).



**For general laboratory use.** Not intended for human or animal diagnostic or therapeutic uses.

**Shipping:** shipped on gel packs

**Storage Conditions:** Store at -20 °C

**Shelf Life:** 12 months