



pLEXY-sat2.1


Integrative constitutive expression vector

Antibiotic selection of transfectants with Nourseothricin (NTC, clonNAT)

Cat. No.	Amount
EGE-274	5 µg

Description:

General expression vector pLEXY-sat2.1 with *sat* marker gene (encoding streptothricine acetyltransferase) allowing selection of recombinant LEXSY strains with the antibiotic LEXSY NTC (Cat.-No. AB-101 and Cat.-No. AB-102).



For general laboratory use. Not intended for human or animal diagnostic or therapeutic uses.

Shipping: shipped on gel packs

Storage Conditions: Store at -20 °C

Shelf Life: 12 months