


**<sup>1,2,8-<sup>13</sup>C<sub>3</sub>-rA-CEP</sup>**[<sup>1,2,8-<sup>13</sup>C</sup>]-N6-Acetyl-5'-O-DMT-2'-O-TBDMS-adenosine-3'-O-CE phosphoramidite

Cat. No.	Amount
N-IT-001-S	10 mg
N-IT-001-L	50 mg

**Applications:**

RNA oligonucleotide synthesis (bottle configured for ABI)

**Isotopic enrichment:** ≥ 98% (Isotopic Enrichment of Precursor)**For general laboratory use.****Shipping:** shipped at ambient temperature**Storage Conditions:** store at -20 °C**Shelf Life:** 12 months after date of delivery**Molecular Formula:** C<sub>45</sub><sup>13</sup>C<sub>3</sub>H<sub>64</sub>N<sub>7</sub>O<sub>8</sub>PSi**Molecular Weight:** 929.12 g/mol**Exact Mass:** 928.44 g/mol**CAS#:** 1220889-84-8 (unlabeled)**Purity:** ≥ 95 % (31P-NMR)**Form:** solid**Color:** white of off-white