


**1,3-¹⁵N₂-rU-CEP**[1,3-¹⁵N]-5'-O-DMT-2'-O-TBDMS-uridine-3'-O-CE phosphoramidite

Cat. No.	Amount
N-IT-007-S	10 mg
N-IT-007-L	50 mg

Applications:

RNA oligonucleotide synthesis (bottle configured for ABI)

Isotopic enrichment: ≥ 98% (Isotopic Enrichment of Precursor)**For general laboratory use.****Shipping:** shipped at ambient temperature**Storage Conditions:** store at -20 °C**Shelf Life:** 12 months after date of delivery**Molecular Formula:** C₄₅H₆₁N₂¹⁵N₂O₉PSi**Molecular Weight:** 863.05 g/mol**Exact Mass:** 862.39 g/mol**CAS#:** 118362-03-1 (unlabeled)**Purity:** ≥ 95 % (31P-NMR)**Form:** solid**Color:** white of off-white