
DATA SHEET





■ NPE-caged-ATP

Adenosine-5'-triphosphate, P3-(1-(2-nitrophenyl)-ethyl)-ester, Sodium salt

Cat. No.	Amount
NU-301S	150 µl (10 mM)
NU-301L	5 x 150 μl (10 mM)

Structural formula of NPE-caged-ATP

For general laboratory use.

Shipping: shipped on gel packs Storage Conditions: store at -20 °C

Additional Storage Conditions: store dark

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

Shelf Life: 12 months after date of delivery Molecular Formula: C₁₈H₂₃N₆O₁₅P₃ (free acid) Molecular Weight: 656.33 g/mol (free acid)

CAS#: 67030-27-7 **Purity:** ≥ 95 % (HPLC) Form: solution in water

Color: colorless to slightly yellow Concentration: 10 mM - 11 mM

Exact Mass: 656.04 g/mol (free acid)

pH: 7.5 ±0.5

Spectroscopic Properties: λ_{max} 260 nm, ϵ 18.0 L mmol⁻¹ cm⁻¹ (Tris-HCl

pH 7.5)

Applications:

Determination of quantum yield and irradiation conditions^[1, 2, 3]

Time resolved solid state NMR^[4]

Ligand for purinergic receptors:

The nucleotide can be transported extra- or intracellular in a protected form to the target. After activation by well-defined conditions the liberated ATP can interact with P2X- and P2Y-receptors. Interacting subreceptor types and corresponding references are listed in Data sheet #NU-1010.

Agonistic ligand, mainly for nucleoside receptor A₁

Nucleoside-triphosphates can be converted by different membranebound phosphatases into nucleosides acting as nucleoside receptor ligands. The caged form is protected during uptake and transport in animal experiments and can be well-directed released through activation at the target tissue.

Selected References:

[1] Omran et al. (2009) Short-length dimethoxynitrophenyl photo-cleavable crosslinkers, synthesis and photolysis. Journal of Photochemistry and Photobiology A: Chemistry 208:125.

[2] Specht et al. (2006) New Photoremovable Protecting Groups for Carboxylic Acids with High Photolytic Efficiencies at Near-UV Irradiation. Application to the Photocontrolled Release of I-Glutamate. ChemBioChem 7:1690.

[3] Zscherp et al. (2001) Reaction-induced infrared difference spectroscopy for the study of protein reaction mechanisms. Biochemistry 40:1875.

[4] Cherepanov et al. (2008) A view on phosphate ester photochemistry by time-resolved solid state NMR. Intramolecular redox reaction of caged ATP. Physical Chemistry Chemical Physics 10:6820.

Volonte et al. (2009) Membrane components and purinergic signalling: the purinome, a complex interplay among ligands, degrading enzymes, receptors and transporters. FEBS J. 276:318.

Yegutkin (2008) Nucleotide and nucleoside converting enzymes: Important modulators of purinergic signalling cascade. Biochim. Biophys. Acta 1783:673.

Scheirlinckx et al. (2001) Monitoring of secondary and tertiary structure changes in the gastric H+/K+-ATPase by infrared spectroscopy Eur. J. Biochem. 268 (13):3644.


Barth et al. (2000) Substrate binding and enzyme function investigated by infrared spectroscopy. FEBS Lett. 477:151.

Hess (1999) Light-Activated (Caged) Biological Ligands. Encyclopedia of Molecular Biology, Vol. 3, T.E. Creighton (Ed.), pp. 1385-1391.

Broustovetsky et al. (1997) Biochemical and physical parameters of the electrical currents measured with the ADP/ATP carrier by photolysis of caged ADP and ATP. Biochemistry 36:13865.



DATA SHEET





■ NPE-caged-ATP

Adenosine-5'-triphosphate, P³-(1-(2-nitrophenyl)-ethyl)-ester, Sodium salt

Higuchi et al. (1997) Kinetics of force generation by single kinesin molecules activated by laser photolysis of caged ATP. Proc. Natl. Acad. Sci. USA 94:4395.

Emoto *et al.* (1995) Tension relaxation induced by pulse photolysis of caged ATP in partially crosslinked fibers from rabbit psoas muscle. *Proc. Natl. Acad. Sci. USA* **92**:1461.

Barth et al. (1995) Photochemical Release of ATP from Caged ATP Studied by Time-Resolved Infrared Spectroscopy. J. Am. Chem. Soc. 117:10311.

Hyman *et al.* (1992) Microtubule-motor activity of a yeast centromere-binding protein complex. *Nature* **359**:533.

Fajer et al. (1990) Myosin heads have a broad orientational distribution during isometric muscle contraction: time-resolved EPR studies using caged ATP. Proc. Natl. Acad. Sci. USA 87:5538.

Walker et al. (1988) Photolabile 1- (2-Nitrophenyl)ethyl Phosphate Esters of Adenine Nucleotide Analogues. Synthesis and Mechanism of

Photolysis. J. Am. Chem. Soc. 110:7170.

Goldman *et al.* (1982) Relaxation of muscle-fibers by photolysis of caged ATP. *Nature* **300**:701.