



N⁶-(6-Aminoheptyl)-ATP-ATTO-647N

N⁶-(6-Aminoheptyl)-adenosine-5'-triphosphate, labeled with ATTO 647N, Triethylammonium salt

Cat. No.	Amount
NU-805-647N	20 µl (1 mM)

For general laboratory use.

Shipping: shipped on gel packs

Storage Conditions: store at -20 °C

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

Shelf Life: 12 months after date of delivery

Molecular Formula: C₅₈H₇₈N₉O₁₅P₃ (free acid)

Molecular Weight: 1234.21 g/mol (free acid)

Exact Mass: 1233.48 g/mol (free acid)

Purity: ≥ 95 % (HPLC)


Form: solution in water

Color: blue


Concentration: 1.0 mM - 1.1 mM

pH: 7.5 ± 0.5

Spectroscopic Properties: λ_{exc} 646 nm, λ_{em} 664 nm, ε 150.0 L mmol⁻¹ cm⁻¹ (Tris-HCl pH 7.5)



Structural formula of N⁶-(6-Aminoheptyl)-ATP-ATTO-647N



excitation and emission spectrum of ATTO 647N