



## 5-Propargylamino-dCTP-DY-776

5-Propargylamino-2'-deoxycytidine-5'-triphosphate, labeled with DY 776, Triethylammonium salt

Cat. No.	Amount
NU-809-776	40 µl (1 mM)

**For general laboratory use.**

**Shipping:** shipped on gel packs

**Storage Conditions:** store at -20 °C

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

**Shelf Life:** 12 months after date of delivery

**Molecular Formula:** C<sub>56</sub>H<sub>65</sub>N<sub>6</sub>O<sub>21</sub>P<sub>3</sub>S<sub>2</sub> (free acid)

**Molecular Weight:** 1315.20 g/mol (free acid)

**Purity:** ≥ 95 % (HPLC)


**Form:** solution in 10 mM Tris-HCl

**Color:** dark blue


**Concentration:** 1.0 mM - 1.1 mM

**pH:** 7.5 ± 0.5

**Spectroscopic Properties:** λ<sub>exc</sub> 771 nm, λ<sub>em</sub> 793 nm, ε 120.0 L mmol<sup>-1</sup> cm<sup>-1</sup> (Tris-HCl pH 7.5)



Structural formula of 5-Propargylamino-dCTP-DY-776



excitation and emission spectrum of DY 776