





## Aminoallyl-UTP-ATTO-488

5-(3-Aminoallyl)-uridine-5'-triphosphate, labeled with ATTO 488, Triethylammonium salt

| Cat. No.   | Amount       |
|------------|--------------|
| NU-821-488 | 20 µl (1 mM) |



excitation and emission spectrum of ATTO 488

**For general laboratory use.**

**Shipping:** shipped on gel packs

**Storage Conditions:** store at -20 °C

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

**Shelf Life:** 12 months after date of delivery

**Molecular Formula:** C<sub>37</sub>H<sub>41</sub>N<sub>6</sub>O<sub>24</sub>P<sub>3</sub>S<sub>2</sub> (free acid)

**Molecular Weight:** 1110.79 g/mol (free acid)

**Exact Mass:** 1110.08 g/mol (free acid)

**Purity:** ≥ 95 % (HPLC)

**Form:** solution in water

**Color:** yellow

**Concentration:** 1.0 mM - 1.1 mM

**pH:** 7.5 ± 0.5

**Spectroscopic Properties:** λ<sub>exc</sub> 500 nm, λ<sub>em</sub> 520 nm, ε 90.0 L mmol<sup>-1</sup> cm<sup>-1</sup> (Tris-HCl pH 7.5)

**Related Products:**

HighYield T7 RNA Synthesis Kit, #RNT-101

HighYield T7 Atto488 RNA Labeling Kit, #RNT-101-488