



5-Propargylamino-CTP-ATTO-665

5-Propargylamino-cytidine-5'-triphosphate, labeled with ATTO 665, Triethylammonium salt

Cat. No.	Amount
NU-831-665	400 µl (0,1 mM)

For general laboratory use.

Shipping: shipped on gel packs

Storage Conditions: store at -20 °C

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

Shelf Life: 12 months after date of delivery

Molecular Formula: C₁₂H₁₉N₄O₁₄P₃ - ATTO-665 (free acid)

Molecular Weight: 1140.22 g/mol (free acid)

Purity: ≥ 95 % (HPLC)


Form: solution in water

Color: blue


Concentration: 0.10 mM - 0.11 mM

pH: 7.5 ± 0.5

Spectroscopic Properties: λ_{exc} 662 nm, λ_{em} 680 nm, ε 160.0 L mmol⁻¹ cm⁻¹ (Tris-HCl pH 7.5)



Structural formula of 5-Propargylamino-CTP-ATTO-665



excitation and emission spectrum of ATTO 665