



### 5-Propargylamino-ddCTP-ATTO-540Q

5-Propargylamino-2',3'-dideoxycytidine-5'-triphosphate, labeled with ATTO 540Q, Triethylammonium salt

Cat. No.	Amount
NU-850-540Q	80 µl (0,5 mM)

**For general laboratory use.**

**Shipping:** shipped on gel packs

**Storage Conditions:** store at -20 °C

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

**Shelf Life:** 12 months after date of delivery

**Molecular Formula:** C<sub>12</sub>H<sub>19</sub>N<sub>4</sub>O<sub>12</sub>P<sub>3</sub> - ATTO 540Q (free acid)

**Molecular Weight:** 1044.22 g/mol (free acid)


**Purity:** ≥ 95 % (HPLC)

**Form:** solution in water


**Concentration:** 0.50 mM - 0.55 mM

**pH:** 7.5 ±0.5

**Spectroscopic Properties:** λ<sub>exc</sub> 543 nm, ε 105.0 L mmol<sup>-1</sup> cm<sup>-1</sup> (Tris-HCl pH 7.5)



Structural formula of 5-Propargylamino-ddCTP-ATTO-540Q



absorption spectrum of ATTO 540Q