



## 5-Propargylamino-ddCTP-ATTO-590

5-Propargylamino-2',3'-dideoxycytidine-5'-triphosphate, labeled with ATTO 590, Triethylammonium salt

Cat. No.	Amount
NU-850-590	160 µl (0,5 mM)

**For general laboratory use.**

**Shipping:** shipped on gel packs

**Storage Conditions:** store at -20 °C

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

**Shelf Life:** 12 months after date of delivery

**Molecular Formula:** C<sub>49</sub>H<sub>55</sub>N<sub>6</sub>O<sub>16</sub>P<sub>3</sub> (free acid)

**Molecular Weight:** 1076.92 g/mol (free acid)

**Purity:** ≥ 95 % (HPLC)


**Form:** solution in water

**Color:** red-violet


**Concentration:** 0.50 mM - 0.55 mM

**pH:** 7.5 ± 0.5

**Spectroscopic Properties:** λ<sub>exc</sub> 593 nm, λ<sub>em</sub> 622 nm, ε 120.0 L mmol<sup>-1</sup> cm<sup>-1</sup> (Tris-HCl pH 7.5)



Structural formula of 5-Propargylamino-ddCTP-ATTO-590



excitation and emission spectrum of ATTO 590