



5-Propargylamino-ddCTP-ATTO-647N

5-Propargylamino-2',3'-dideoxycytidine-5'-triphosphate, labeled with ATTO 647N, Triethylammonium salt

Cat. No.	Amount
NU-850-647N	80 µl (0,5 mM)

For general laboratory use.

Shipping: shipped on gel packs

Storage Conditions: store at -20 °C

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

Shelf Life: 12 months after date of delivery

Molecular Formula: C₅₄H₆₈N₇O₁₄P₃ (free acid)

Molecular Weight: 1132.08 g/mol (free acid)

Exact Mass: 1131.40 g/mol (free acid)

Purity: ≥ 95 % (HPLC)


Form: solution in water

Color: blue


Concentration: 0.50 mM - 0.55 mM

pH: 7.5 ± 0.5

Spectroscopic Properties: λ_{exc} 646 nm, λ_{em} 664 nm, ε 150.0 L mmol⁻¹ cm⁻¹ (Tris-HCl pH 7.5)



Structural formula of 5-Propargylamino-ddCTP-ATTO-647N



excitation and emission spectrum of ATTO 647N