



8-(6-Aminoethyl)-amino-cAMP-ATTO-633

8-(6-Aminoethyl)-amino-adenosine-3',5'-cyclic monophosphate, labeled with ATTO 633, Triethylammonium salt

Cat. No.	Amount
NU-851-633	80 µl (1 mM)

For general laboratory use.

Shipping: shipped on gel packs

Storage Conditions: store at -20 °C

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

Shelf Life: 12 months after date of delivery

Molecular Formula: C₅₁H₆₅N₁₀O₈P (free acid)

Molecular Weight: 977.12 g/mol (free acid)

Exact Mass: 976.47 g/mol (free acid)

Purity: ≥ 95 % (HPLC)


Form: solution in water

Color: blue


Concentration: 1.0 mM - 1.1 mM

pH: 7.5 ± 0.5

Spectroscopic Properties: λ_{exc} 630 nm, λ_{em} 651 nm, ε 130.0 L mmol⁻¹ cm⁻¹ (Tris-HCl pH 7.5)



Structural formula of 8-(6-Aminoethyl)-amino-cAMP-ATTO-633



excitation and emission spectrum of ATTO 633