


**cyclic-di-AMP**

(c-di-AMP)

Cyclic diadenosine monophosphate, Sodium salt

3',5'-Cyclic diadenylic acid sodium salt

Cat. No.	Amount
NU-954S	1 μ mol
NU-954L	5 x 1 μ mol



Structural formula of cyclic-di-AMP

For general laboratory use.**Shipping:** shipped on gel packs**Storage Conditions:** store at -20 °C

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

Shelf Life: 12 months after date of delivery**Molecular Formula:** C₂₀H₂₄N₁₀O₁₂P (free acid)**Molecular Weight:** 658.41 g/mol (free acid)**Exact Mass:** 658.11 g/mol (free acid)**CAS#:** 54447-84-6**Purity:** \geq 98 % (HPLC)**Form:** solid**Color:** white to off-white**Spectroscopic Properties:** λ_{max} 259 nm, ϵ 27.0 L mmol⁻¹ cm⁻¹ (Tris-HCl pH 7.5)