

**HEV ORF2-ORF3 32 kDa**

Hepatitis E virus ORF2 -ORF3 protein, mosaic recombinant, *E. coli*

| Cat. No.  | Amount |
|-----------|--------|
| PR-1420-S | 200 µg |
| PR-1420-L | 1 mg   |

**For general laboratory use.**

**Shipping:** shipped on gel packs

**Storage Conditions:** store at -20 °C

**Additional Storage Conditions:** avoid freeze/thaw cycles

**Shelf Life:** 6 months

**Molecular Weight:** about 32 kDa

**Form:** liquid (Supplied in Tris-HCl 50 mM pH 8.0, Sodium Chloride 60 mM, Glutathione 10 mM, Urea 0.5 M, Glycerol 50 %)

**Applications:**

Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Description:**

The protein contains immunodominant regions of HEV ORF2 and ORF3, amino acids (81-133)(631-660). The protein is fused to a GST tag.

**Specificity:**

Immunoreactive with sera of HEV-infected individuals