

**Legionella pneumophila Antigen**

Pool of Legionella pneumophila serovars 1A, 1B, 2, 3, 4, 5, 6 and 7 native, Solid medium

Cat. No.	Amount
PR-BA106-S	250 µg
PR-BA106-L	1 mg

For general laboratory use.

Shipping: shipped on dry ice

Storage Conditions: store at -80 °C

Additional Storage Conditions: avoid freeze/thaw cycles

Shelf Life: 12 months

Form: liquid (Supplied in PBS)

Culture system: Solid medium

Protein concentration: Lot specific. Determined by Bradford assay (Biorad) using BSA as standard.

Preservation: None

Analysis: ELISA

Preparation: Washed and heat treated organisms are presented as a suspension in phosphate buffered saline.

Inactivation: This product has been inactivated by heat treatment. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.