

**Parvovirus B19 VP2 VLPs**

recombinant, Sf9 insect cells

Cat. No.	Amount
PR-BA122VP2-S	250 µg
PR-BA122VP2-L	1 mg

For general laboratory use.**Shipping:** shipped on dry ice**Storage Conditions:** store at -80 °C**Additional Storage Conditions:** avoid freeze/thaw cycles**Shelf Life:** 12 months**Form:** liquid**Concentration:** Lot specific, specified in mg/ml. Determined by Bradford assay (Biorad) using BSA as standard.**Description:**

Virus-like particles from VP2 structural proteins of parvovirus B19.

Purification:

Chromatography

Presentation:

The antigen is presented in phosphate buffered saline at pH 7.4, containing 10 % glycerol.

Analysis:

SDS-PAGE and ELISA