

***Yersinia enterocolitica* O:8 Antigen**

native, liquid culture

Cat. No.	Amount
PR-BA13808-S	250 µg
PR-BA13808-L	1 mg

For general laboratory use.**Shipping:** shipped on dry ice**Storage Conditions:** store at -80 °C**Additional Storage Conditions:** avoid freeze/thaw cycles**Shelf Life:** 12 months**Form:** liquid (Supplied in PBS)**Culture system:** Liquid culture**Protein concentration:** Lot specific. Determined by Bradford assay (Biorad) using BSA as standard.**Preservation:** None**Analysis:** SDS-PAGE and ELISA**Preparation:** *Yersinia* is grown in liquid culture and stimulated to secrete YOPs (*Yersinia* outer proteins). YOP proteins concentrated from the culture medium are presented in phosphate buffered saline.**Inactivation:** This product has been inactivated by sterile filtration. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.