

**Quantifoil™ R 0.6/1 on 300 copper mesh**

Cat. No.	Amount
X-100-CU300	100 pcs.

For general laboratory use.

Shipping: shipped at ambient temperature

Storage Conditions: store at ambient temperature

Shelf Life: n/a

Applications:

Support film for biological samples in Cryo-EM techniques

Description:

QUANTIFOIL™ R 0.6/1 carbon film has circular holes of 0.6 µm diameter with an interspace of 1 µm in both dimensions, applied on 300 copper mesh.