

**Quantifoil™ R 2/1 on 300 gold mesh**

Cat. No.	Amount
X-102-Au300	100 pcs.

For general laboratory use.

Shipping: shipped at ambient temperature

Storage Conditions: store at ambient temperature

Shelf Life: n/a

Applications:

Support film for biological samples in Cryo-EM techniques

Description:

QUANTIFOIL™ R 2/1 carbon film has circular holes of 2 µm diameter with an interspace of 1 µm in both dimensions, applied on biocompatible 300 gold mesh. Suitable for magnifications between 30.000x and 40.000x.