



## C-flat™ 1.2/1.3 on 200 gold mesh

| Cat. No.    | Amount  |
|-------------|---------|
| X-301-AU200 | 50 pcs. |


### Applications:

Support film for biological samples in Cryo-EM techniques

### Description:

20 nm C-flat™ carbon film with 1.2 µm hole size and 1.3 µm hole spacing, applied on 200 gold mesh.

**nonsaleable to USA**



**For general laboratory use.**

**Shipping:** shipped at ambient temperature

**Storage Conditions:** store at ambient temperature

**Shelf Life:** n/a