



## C-flat™ 1.2/1.3 on 300 gold mesh

Cat. No.	Amount
X-301-AU300	50 pcs.

**For general laboratory use.**

**Shipping:** shipped at ambient temperature

**Storage Conditions:** store at ambient temperature

**Shelf Life:** n/a

**Applications:**

Support film for biological samples in Cryo-EM techniques

**Description:**

20 nm C-flat™ carbon film with 1.2 µm hole size and 1.3 µm hole spacing, applied on 300 gold mesh.

**nonsaleable to USA**