



**C-flat™ 1.2/1.3 on 200 copper mesh**

Cat. No.	Amount
X-301-CU200	50 pcs.


**Applications:**

Support film for biological samples in Cryo-EM techniques

**Description:**

20 nm C-flat™ carbon film with 1.2 µm hole size and 1.3 µm hole spacing, applied on 200 copper mesh.

**nonsaleable to USA**



**For general laboratory use.**

**Shipping:** shipped at ambient temperature

**Storage Conditions:** store at ambient temperature

**Shelf Life:** n/a