





## C-flat™ 2/1 on 200 gold mesh

Cat. No.	Amount
X-302-AU200	50 pcs.



**Carbon Film**  
20 nm thick  
2.0 µm holes  
1.0 µm spacing

**Gold Grid**  
200 Mesh

**Pack Size:** 50 Grids

**Applications:**

Support film for biological samples in Cryo-EM techniques

**Description:**

20 nm C-flat™ carbon film with 2 µm hole size and 1 µm hole spacing, applied on 200 gold mesh.

**nonsaleable to USA**

**For general laboratory use.**

**Shipping:** shipped at ambient temperature

**Storage Conditions:** store at ambient temperature

**Shelf Life:** n/a