

**C-flat™ 2/1 on 300 copper mesh**

| Cat. No. | Amount |
|-------------|---------|
| X-302-CU300 | 50 pcs. |

For general laboratory use.

Shipping: shipped at ambient temperature

Storage Conditions: store at ambient temperature

Shelf Life: n/a

Applications:

Support film for biological samples in Cryo-EM techniques

Description:

20 nm C-flat™ carbon film with 2 µm hole size and 1 µm hole spacing, applied on 300 copper mesh.

nonsaleable to USA