

**C-flat™ 2/2 on 400 gold mesh**

| Cat. No. | Amount |
|-------------|---------|
| X-303-AU400 | 50 pcs. |

For general laboratory use.

Shipping: shipped at ambient temperature

Storage Conditions: store at ambient temperature

Shelf Life: n/a

Applications:

Support film for biological samples in Cryo-EM techniques

Description:

20 nm C-flat™ carbon film with 2 µm hole size and 2 µm hole spacing, applied on 400 gold mesh.

nonsaleable to USA