

**Au-flat™ 0.6/1.0 on 300 gold mesh**

Cat. No.	Amount
X-400-AU300	50 pcs.

**For general laboratory use.**

**Shipping:** shipped at ambient temperature

**Storage Conditions:** store at ambient temperature

**Shelf Life:** n/a

**Applications:**

Support film for biological samples in Cryo-EM techniques

**Description:**

45 nm gold alloy film (80% Au and 20% Pd) with 0.6 µm hole size and 1.0 µm hole spacing, applied on 300 gold mesh

**nonsaleable to USA**