



## Cryo-EM Grid Boxes

Cat. No.	Amount
X-CEM-201	20 pcs.



Swissci Cryo-EM grid box (label not included)



Cryo-EM puck with 11 labeled Swissci Cryo-EM grid boxes

**For general laboratory use.**

**Shipping:** shipped at ambient temperature

**Storage Conditions:** store at ambient temperature

**Shelf Life:** n/a

**Description:**

Compatible with common Cryo-EM sample mounting and storage devices, the Swissci Cryo-EM grid box offers secure Cryo-EM sample storage.

The round Cryo-EM grid box, made from an anti-static polymer, has four diamond-shaped slot positions that each hold one Cryo-EM grid, e.g. QUANTIFOIL™ Holey Carbon Film. The grid box can be labeled with any cryo label/marker. The Cryo-EM samples can be identified by the unique grid box label and the slot position (I-III) indicated on the side of the box. Up to 12 Swissci Cryo-EM grid boxes can be stored together in a Cryo-EM Puck.

**Related Products:**

Cryo-EM Puck Gen 2.0, #X-CEMG2

Flat Tip Tweezers, #X-TW-2